PURACOL®
PURACOL® PLUS
INSTRUCTIONS FOR USE

Puracol® or Puracol® Plus MircroScaffold™ wound dressings are advanced wound dressings composed of pure collagen. This unique, sterile biomaterial forms a soft, conformable moist gel sheet at the wound surface, promoting natural healing. The collagen based, moist microenvironment provides an ideal scaffold for the infiltration, proliferation and growth of cells involved in natural wound healing, allowing optimal levels of granular tissue formation and rapid epithelialization. The porous microstructure of the product promotes capillary action, which in turn adds to the intrinsic absorbency of native collagen. The absorbent nature of the dressing further allows for the removal of necrotic tissue fragments, exudate, and miscellaneous debris from the wound site. The dressing can be used as a primary dressing, or be used in combination with other traditional or advanced wound dressings of occlusive or semi-occlusive nature. The dressing can be cut to any size to fit any wound.

INDICATIONS:
- Full thickness and partial thickness wounds
- Pressure and venous ulcers
- Ulcers caused by mixed vascular etiologies
- Diabetic ulcers
- Partial and full-thickness burns
- Donor sites and other bleeding surface wounds
- Abrasions
- Traumatic wound healing by secondary intention
- Dehisced surgical incisions

GENERAL PRECAUTIONS:
Puracol® or Puracol® Plus dressings are suitable for use under compression under proper supervision by a healthcare professional. Wounds with infection may be managed with Puracol® or Puracol® Plus dressings when the causes of infection are being treated in parallel by medical professionals.

CONTRAINDICATIONS:
Puracol® or Puracol® Plus cannot be used on wounds with active vasculitis or patients with known sensitivity to collagen; product not indicated for third-degree burns.

DIRECTIONS FOR USE:
1. Follow normal procedures for wound bed preparation such as irrigation with wound cleansers or sterile saline, debriding the wound if necessary.
2. Dry the area surrounding the wound.
3. The product should be directly applied to the wound, and may be moistened with sterile saline to pre-form the gel if the wound bed is very dry.
4. The dressing may be cut to size or shape to better match the wound area.
5. The product should be covered with a suitable secondary dressing in order to maintain the moist wound environment.
6. The dressing may be replaced at the discretion of the healthcare professional. The product is 100% pure collagen and is completely biodegradable. During dressing changes, it is not necessary to remove any residual Puracol or Puracol Plus.

SHELF LIFE:
When stored properly, Puracol® or Puracol® Plus collagen wound dressing has a shelf life of 3 years. If Puracol® or Puracol® Plus is not applied immediately after opening the primary pouch, it should be discarded. Do not attempt to re-sterilize the dressing. Do not use the product after its expiration date has passed.

STORAGE:
Store in a cool, dry place. Avoid excessive heat and humidity.

Single use only. Rx only.
Sterility guaranteed in unopened, undamaged package.
This product is not made with Natural Rubber Latex.

MSC8511	Puracol, 1 x 1 in (2.5 x 2.5 cm)	10/bx, 5 bx/cs
MSC8522BB	Puracol, 2 x 2 in (5.1 x 5.1 cm)	30/bg, 10 bg/cs
MSC8522	Puracol, 2 x 2 in (5.1 x 5.1 cm)	10/bx, 5 bx/cs
MSC8544	Puracol, 4 x 4.25 in (10.1 x 10.8 cm)	10/bx, 5 bx/cs
MSC8588	Puracol, 8 x 8 in (20.3 x 20.3 cm)	10/bx, 5 bx/cs
MSC861X8EP	Puracol Plus, 1 x 8 in (2.5 x 20 cm)	10/bx, 5 bx/cs
MSC8622EP	Puracol Plus, 2 x 2.2 in (5 x 5.6 cm)	10/bx, 5 bx/cs
MSC8644EP	Puracol Plus, 4.2 x 4.5 in (10.7 x 11.4 cm)	10/bx, 5 bx/cs

www.medline.com ©2020 Medline and Puracol are registered trademarks and MicroScaffold is a trademark of Medline Industries, LP. Made in Germany by MedSkin Solutions Dr. Suwelack for Medline Industries, LP, Three Lakes Drive, Northfield, IL 60093 USA. **1-800-MEDLINE** V4 RE22DSS

medskin solutions DR. SUWELACK