

EFFECTIVE DATE: 10/29/86 REVISION OF: 11/15/02

# DISPATCH<sup>®</sup> Material Safety Data Sheet

## 1. Product Identity

### DISPATCH<sup>®</sup> Hospital Cleaner Disinfectant with Bleach

#### ITEM NUMBERS:

68963 — 8 oz Pump Spray Bottle  
 68967 — 22 oz Trigger Spray Bottle  
 68970 — 32 oz Trigger Spray Bottle  
 68832 — 32 oz Pull Top Bottle  
 68973 — 64 oz Refill  
 68978 — 128 oz Refill  
 69260 — 60 Piece Soft Pack\*  
 69101 — 50 Single-Piece Towels\*

\*Safety information pertains to solution on towels.

Manufactured by: Caltech Industries, Inc.  
 2420 Schuette Drive, Midland, MI 48642

## 2. Ingredients

(Typical Values — Not Specifications)

Active Ingredients	CAS NO.	%	OSHA PEL	ACGIH TLV
Sodium Hypochlorite	007681-52-9	<1	Not Est.	Not Est.
Inert Ingredients <sup>1</sup>	--	>99	—	—

<sup>1</sup>Includes water, detergents, stabilizer and fragrance

## 3. Physical Data

**BOILING POINT:** 212°F, 100°C.

**VAP PRESSURE:** 17.5 mmHg at 20°C.

**VAP DENSITY:** Not determined.

**APPEARANCE:** Clear, colorless solution.

**SOLUBILITY IN WATER:** Completely soluble.

**SPECIFIC GRAVITY:** 1.013-1.017 (25°/25°C).

**pH:** 12.2.

**ODOR:** Slight bleach odor.

## 4. Fire & Explosion Hazard Data

**NFPA HAZARD RATING:** 0

**FLASH POINT:** None

**FLAMMABLE LIMITS—LFL/UFL:** Not applicable

**METHOD USED:** TCC

**EXTINGUISHING MEDIA:** Not combustible

**SPECIAL FIRE FIGHTING EQUIPMENT & HAZARDS:** None.

## 5. Reactivity Data

**NFPA HAZARD RATING: 1**

**STABILITY (Conditions To Avoid):** Stability decreases with heat, light exposure, decreases in pH and contamination with heavy metals, such as nickel, cobalt, copper and iron.

**INCOMPATIBILITY (Specific Materials To Avoid):** Avoid contact with acid, reducing agents, or ammonia.

**HAZARDOUS DECOMPOSITION PRODUCTS:** Chlorine in the presence of acids.

**HAZARDOUS POLYMERIZATION:** Will not occur.

**6. Specific Hazards****NFPA HAZARD RATING: OX**

Material is basic and an oxidizer. Do not mix with acidic materials such as certain toilet cleaners, ammonia, or with formaldehyde.

**7. Health Hazard Data****NFPA HAZARD RATING: 0**

**EYES:** May cause moderate irritation with corneal injury.

**SKIN CONTACT:** Prolonged or repeated exposure not likely to cause significant skin irritation.

**SKIN ABSORPTION:** A single prolonged exposure is not likely to result in the material being absorbed through the skin in harmful amounts. The LD50 for skin absorption in rabbits is >4 g/kg.

**INGESTION:** Single dose oral toxicity is low. The oral LD50 for male and female rats is >5 g/kg. Amounts ingested incidental to industrial handling are not likely to cause injury.

**INHALATION:** No adverse effects are anticipated from inhalation.

**SYSTEMIC & OTHER EFFECTS:** No specific data available. However, repeated exposures are not anticipated to cause any significant adverse effects.

**CARCINOGENICITY:** Not listed as carcinogen by NTP, OSHA or IARC.

**8. First Aid**

**EYES:** Irrigate with flowing water immediately and continuously for 15 minutes. Consult medical personnel.

**SKIN:** Wash off in flowing water or shower.

**INGESTION:** Drink a glass of water. Consult medical personnel.

**INHALATION:** No adverse effects anticipated by this route of exposure incidental to proper handling.

**NOTE TO PHYSICIAN:** No specific antidote. Supportive care. Treatment based on judgment of the physician in response to reactions of the patient.

**9. Special Handling Information**

**VENTILATION:** Good general ventilation should be sufficient for most conditions.

**RESPIRATORY PROTECTION:** No respiratory protection should be needed.

**PROTECTIVE CLOTHING:** No precautions other than clean body-covering clothing should be needed.

**10. Spill, Leak & Disposal Procedures**

**ACTION TO TAKE FOR SPILLS/LEAKS:** Isolate damaged goods. Spills should be diluted with water, then soaked up with suitable absorbent material.


**DISPOSAL METHOD:** Landfill in accordance with local, state, and federal regulations.

**11. Special Precautions & Additional Information**

**PRECAUTIONS TO BE TAKEN IN HANDLING AND STORAGE:** Do not mix with acidic materials such as certain toilet cleaners, ammonia, or with formaldehyde. Avoid eye contact.

The information herein is given in good faith, but no warranty, expressed or implied, is made.

DISPATCH is a registered trademark of Caltech Industries, Inc.



**EMERGENCY PHONE: (989) 496-3110**

## Product Information

# DISPATCH<sup>®</sup> Hospital Cleaner Disinfectant with Bleach

## Rapid-Acting, One-Step Cleaner Disinfectant

DISPATCH<sup>®</sup> Hospital Cleaner Disinfectant with Bleach is a rapid-acting cleaner disinfectant that works in one step. A unique, stabilized bleach and detergent solution (bleach solution strength similar to the 1:10 solution strength recommended by the CDC for effective disinfecting), DISPATCH<sup>®</sup> offers bleach strength disinfecting without the time-consuming task of mixing bleach; it also offers superior cleaning capability.

## Fast-Acting Efficacy, EPA-Registered

This product has demonstrated efficacy for a broad range of organisms in one minute. The results of efficacy tests, conducted in accordance with the test methods of the Association of Official Analytical Chemists in the presence of organic soil (5% blood serum), are described below:

Organism	Contact Time
<i>Mycobacterium tuberculosis</i> (TB) (Two minutes for DISPATCH <sup>®</sup> Towels)	30 seconds
<i>Staphylococcus aureus</i>	1 minute
<i>Salmonella choleraesuis</i>	1 minute
<i>Streptococcus pyogenes</i>	1 minute
<i>Pseudomonas aeruginosa</i>	1 minute
Methicillin resistant <i>Staphylococcus aureus</i> (MRSA) (DISPATCH <sup>®</sup> Liquid only)	1 minute
Vancomycin resistant <i>Enterococcus faecalis</i> (VRE) (DISPATCH <sup>®</sup> Liquid only)	1 minute
HIV-1 (AIDS)	1 minute
Poliovirus (Type-1, Mahoney)	1 minute
Herpes Simplex (HSV-2)	1 minute
Rotavirus (DISPATCH <sup>®</sup> Liquid only)	1 minute
<i>Trichophyton mentagrophytes</i>	1 minute

## Meets Bloodborne Pathogen Guidelines

This product meets OSHA, CDC, NCCLS and Universal Precautions Guidelines for healthcare worker protection from bloodborne pathogen contamination.

## Cleans & Disinfects in One Step

- Excellent for removing blood and blood products, body fluids and excretions and everyday soil

- Cleans and disinfects in one step; no precleaning required; rarely needs rinsing

### **Convenient**

- Available in liquid (sprays or refills) and pre-moistened, disposable towels

### **Versatile, Multiple Use Product**

- Especially useful for laboratory countertops, point-of-care testing equipment such as blood glucose monitors, IV poles, respiratory and physical therapy equipment as well as general surfaces such as carts, bed frames and springs
- Excellent for cleaning and disinfecting laminated plastic countertops, stainless steel, plastic, ceramic tile, grout, vinyl plastics, fiberglass and porcelain

### **Outperforms Conventional Bleach Solutions**

Unlike hand-diluted bleach solutions, this product:

- Is convenient and ready-to-use; no mixing or discarding of unused portions
- Does not rapidly deteriorate; stays stable to expiration date
- Is economical; eliminates daily mixing and use of a separate cleaning product
- Contains a pleasant fragrance and anti-corrosive agent

### **Size & Order Information**

Item No.

68963 — 8 oz spray, (15) per case

68967 — 22 oz spray, (8) per case

68970 — 32 oz spray, (6) per case

68832 — 32 oz pull top, (6) per case

68973 — 64 oz refill, (6) per case

68978 — 128 oz refill, (4) per case

69260 — 60 piece soft pack, (12) per case

69101 — 50 single-piece towels and dispenser, (6) dispensers per case

**EPA Registration No. 56392-7**

**EPA Registration No. 56392-8**



---

**CALTECH**  
**INDUSTRIES, INC.**

---

2420 Schuette Drive  
Midland, MI 48642