

**NEW!**

# microdot® BLEACH WIPE

**Can be used on blood glucose meters.**

## KILL CLAIMS:

- Hepatitis B Virus (HBV)
- Methicillin-Resistant *Staphylococcus aureus* (MSRA)
- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- Vancomycin-Resistant *Enterococcus faecalis* (VRE)
- Human Immunodeficiency Virus Type 1 (HIV-1)
- *Trichophyton mentagrophytes*

## PRODUCT DESCRIPTION:

Microdot's® Bleach (sodium hypochlorite) saturated wipes in a pop-up canister are EPA registered, hospital use disinfectant wipes. These wipes contain 0.525% sodium hypochlorite (equivalent to the "1:10 dilution." This product is designed for cleaning, decontaminating, and disinfecting hard, non-porous healthcare surfaces and equipment.

Keep these wipes in a tightly closed canister when not in use. Store these products in a cool, dry, and well-ventilated area away from heat, open flame, direct sunlight and inaccessible to children. Rotate dated products while in inventory/storage.

## SUGGESTED AREAS OF USAGE:

- Hospitals
- Nursing Homes
- Doctor's Offices
- Dental Offices
- Point-of-Care Equipment
- Emergency Rooms
- Laboratories
- EMS
- Out-Patient Care
- Surgery
- Environmental Services
- Infection Control

**160 Wipes per canister • 6" x 6" wipe  
12 canister/case (1,920 wipes)**

## EPA Registered, Ready to Use hospital use 1:10 dilution

**Can be used to disinfect hard, non-porous surfaces in hospitals, doctors' offices, nursing homes, schools, and dental offices.**

**ASK YOUR DISTRIBUTOR FOR  
Microdot® Bleach Wipe**



**GROVE MEDICAL**

**Local: 864.269.0283**

**Toll FREE: 800.269.1405**

**www.GroveMedical.com**

**ITEM NUMBER: 600-12**



**CAMBRIDGE SENSORS USA, LLC**

**www.microdotcs.com**