


Shiley™ tracheostomy tube accessories



8A. Disposable inner cannula (DIC) with integral 15 mm snap-lock connector, individually sterile packaged, required for use with DCT, DFEN, PERC, DCFS, and DCFN.

8B. Spare inner cannula (SIC) with red, color-coded twist-lock connector, for temporary use during cleaning or removal of the inner cannula of Shiley™ tracheostomy tubes, may be used with LPC, FEN, CFS, and CFN.

8C. White 15 mm cap (CAP) attaches to the fenestrated inner cannula with green 15 mm connector to permit breathing through the fenestration and upper airway.

8D. Decannulation plug (DCP) closes the proximal end of the FEN and CFN tubes to permit breathing through the fenestration and upper airway.

8E. Disposable decannulation plug (DDCP) closes the proximal end of the DCFS, DCFN, and DFEN to permit breathing through the fenestration and upper airway.

SIZE	PRODUCT DESIGNATION	FOR USE WITH			QTY PER BOX	
4	4DIC		4DFEN	4DCFN	10	DIC
6	6DIC		6DFEN	6DCFN	10	
8	8DIC		8DFEN	8DCFN	10	
10	10DIC		10DFEN	10DCFN	10	
4	4SIC		4FEN	4CFS	10	SIC
6	6SIC		6FEN	6CFS	10	
8	8SIC		8FEN	8CFS	10	
10	10SIC		10FEN	10CFS	10	
4	4DICFEN	4DFEN	4DCFN		10	DICFEN
6	6DICFEN	6DFEN	6DCFN		10	
8	8DICFEN	8DFEN	8DCFN		10	
10	10DICFEN	10DFEN	10DCFN		10	
4	4DCP	4FEN	4CFN		10	DCP
6	6DCP	6FEN	6CFN		10	
8	8DCP	8FEN	8CFN		10	
10	10DCP	10FEN	10CFN		10	
	DDCP	Any size DCFS, DCFN, and DFEN tracheostomy tube			10	DDCP
	CAP	Any size FEN and CFN tracheostomy tube			10	CAP