


Rapid Tests Pregnancy

QuPID® hCG

Single step cassette formats



>99% accurate / >99% specific

Room temperature storage


Built-in positive and negative procedural controls

CE marked

QuPID® hCG	Reference No.	Contents
QuPID® hCG 20 mIU/mL in urine	1220-025	25 tests
	1220-050	50 tests
QuPID® Plus hCG 20 mIU/mL in urine or 10 mIU/mL in serum	1230-025	25 tests
	1230-050	50 tests

True®20 hCG

Single step cassette formats



>99% accurate / >99% specific

Room temperature storage

Built-in positive and negative procedural controls

CE marked

True®20 hCG	Reference No.	Contents
True® 20 hCG 20 mIU/mL in urine	1430-050	50 tests
True® 20 Plus hCG 20 mIU/mL in urine or serum	1440-050	50 tests

Pregnancy test controls

CE marked

Controls	Reference No.	Contents
hCG Bi-level urine controls (High positive and negative)	1225-205	2 x 5 mL bottles
hCG Tri-level serum controls (Low positive, high positive and negative)	1235-305	3 x 5 mL bottles