

Tamper-Proof Magnetic Pull Cord Alarm

HCPCS:



- When magnetic pull switch is removed, alarm is activated
- Pull switch can be activated from any direction
- Convenient On/Off switch for easy activation and deactivation
- Easily mounts on wheelchair or bed
- Cord adjusts from 28" 58"
- Alligator clip allows for attachment to patients clothes
- 97-103 dB
- One 9V battery included

PRODUCTS

Description

PRODUCTS	Description
13603	Pull Cord Alarm,w/Magnet,1/ea

ACCESSORIES

Description

ACCESSORIES	Description
136DC	ALARM AC ADAPTER, 1/EA, 25/CV, CS
13617	Cord,Magnet,Clip for 13603,1/ea, EA