

## SenSura® Mio Convex Light MAXI Urostomy Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



**Transparent** 10<sup>1</sup>/<sub>2</sub>" (27 cm) / 400 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
16835	Cut-to-fit	$\frac{3}{8}$ - $\frac{7}{8}$ " (10-23 mm)	A4392
16836	Cut-to-fit	$\frac{5}{8}$ - $1\frac{5}{16}$ " (15-33 mm)	A4392
16837	Cut-to-fit	$\frac{3}{8}$ - $1\frac{11}{16}$ " (10-43 mm)	A4392
16840	Pre-cut	$1\frac{3}{16}$ " (21 mm)	A4392
16841	Pre-cut	1" (25 mm)	A4392
16842	Pre-cut	$1\frac{1}{8}$ " (28 mm)	A4392

**Opaque** 10<sup>1</sup>/<sub>2</sub>" (27 cm) / 400 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
16825	Cut-to-fit	$\frac{3}{8}$ - $\frac{7}{8}$ " (10-23 mm)	A4392
16826	Cut-to-fit	$\frac{5}{8}$ - $1\frac{5}{16}$ " (15-33 mm)	A4392
16830	Pre-cut	$1\frac{1}{16}$ " (18 mm)	A4392
16831	Pre-cut	$1\frac{3}{16}$ " (21 mm)	A4392
16832	Pre-cut	1" (25 mm)	A4392
16833	Pre-cut	$1\frac{1}{8}$ " (28 mm)	A4392
16827	Pre-cut	$1\frac{1}{4}$ " (31 mm)	A4392

## SenSura® Mio Convex Light MIDI Urostomy Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Multi-chambers minimize sloshing
- Anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Neutral gray, water-resistant textile
- **Not made with natural rubber latex**



**Opaque** 9<sup>1</sup>/<sub>2</sub>" (24 cm) / 330 mL Units: 10

Code	Cut	Stoma Size	HCPCS*
16821	Cut-to-fit	$\frac{3}{8}$ - $\frac{7}{8}$ " (10-23 mm)	A4392
16822	Cut-to-fit	$\frac{5}{8}$ - $1\frac{5}{16}$ " (15-33 mm)	A4392
16823	Pre-cut	$1\frac{1}{4}$ " (31 mm)	A4392

 Indicates convexity

\*For important information regarding the use of HCPCS codes listed and code definitions see page 192