RESTING HAND SPLINT



Product Information

Measurement*	Size	Right	Left
21/4" - 23/4"	Small	325BR	325BL
21/8" - 31/4"	Medium	325CR	325CL
33/8" - 33/4"	Large	325DR	325DL

^{*} Measure width of hand at MP joints

Indications

- Post-surgery/injury positioning, wrist drop and finger flexor tightness
- Not recommended when moderate to severe spasticity is present

Features & Benefits

- Supports the fingers, thumb and wrist in functional/resting position
- Malleable aluminum frame allows for static progressive splinting
- Simple adjustment of the wrist and hand with no tools or heat required
- Self-wicking liner helps control perspiration and dampness to prevent skin breakdown
- Hand washable



Powell, TN 37949 www.deroyal.com 888.938.7828 • 865.938.7828 rev 11/07