



THIN



PLUS



PLUS V.A.D.

COVADERM® island dressing



moist wound

- Non-adherent, general purpose sterile dressing
- Protective, highly absorbent sterile pad
- Non-irritating acrylic adhesive tape
- Comfortable rounded edges and stretchable fabric tape
- Split release liner for easy one-step application
- Conforms securely to body contours

COVADERM PLUS® adhesive composite dressing



moist wound

- Multi-layered composite wound dressing
- Helps maintain moist environment with protection from outside contaminants
- Four-layer construction:
 - Protective non-adherent contact layer
 - Soft pad for drainage absorption
 - Breathable, semi-occlusive polyurethane barrier film
 - Conformable adhesive-coated fabric tape
- Helps prevent drainage strike-through
- Helps prevent wound trauma
- Ideal for many wound types
- Available with radial slit for tube sites

Product #	Overall	Pad	Qty/Bx	Bx/Cs	Qty/Cs	HCPC
COVADERM®						
46-001	4" x 4"	2½" x 2½"	25	8	200	A6254
46-002	4" x 6"	2½" x 4"	25	8	200	A6254
46-003	4" x 8"	2" x 5½"	25	4	100	A6254
46-004	4" x 10"	2" x 7½"	25	4	100	A6254
46-005	4" x 14"	2" x 11"	25	4	100	A6255
COVADERM® THIN						
15-201	4" x 4"	2½" x 2½"		100		A6219
15-202	4" x 6"	2½" x 4"		100		A6219
15-203	4" x 8"	2½" x 5½"		50		A6219
15-204	4" x 10"	2½" x 7½"		50		A6219
15-205	4" x 14"	2½" x 11"		50		A6220

indications:

Open and minimally invasive surgical incision sites
Arthroscopic incisions
Superficial lacerations
Abrasions

Chronic wounds
Pressure ulcers
Leg ulcers
Laparotomy sites
Diabetic ulcers

contraindications:

Do not use on:
third-degree burns, heavily bleeding wounds, individuals

who are sensitive to or have had an allergic reaction to the dressing or its materials



Product #	Tape	Pad	Qty/Bx	Qty/Cs	HCPC	
COVADERM PLUS®						
46-400	2" x 2"	1" x 1"	25	100	A6203	
46-411	2½" x 2½"	1½" x 1½"	10	50	A6203	
46-401	4" x 4"	2½" x 2½"	10	50	A6203	
46-402	6" x 6"	4" x 4"	10	50	A6203	
46-403	6" x 8"	4" x 6"	10	50	A6204	
46-414	4" x 10"	2" x 7½"	N/A	100	A6203	
46-415	4" x 14"	2" x11"	N/A	100	A6204	
46-408	4" Diameter w/2" Radial Slit	2½" Diameter	25	200	A6203	
COVADERM PLUS® V.A.D.						
46-405	4" x 4" V.A.D.	2" x 2"	1" x 4"	25	200	A6255
46-406	6" x 6" V.A.D.	4" x 4"	1" x 6"	5	50	A6255
COVADERM PLUS® THIN						
15-301	2" x 2"	1" x 1"	50	100/cs	A6203	
15-302	4" x 4"	2½" x 2½"	50	100/cs	A6203	
15-303	6" x 6"	4" x 4"	50	50/cs	A6203	
15-304	6" x 8"	4" x 6"	50	50/cs	A6204	
15-305	4" x 10"	2" x 7½"	50	50/cs	A6203	
15-306	4" x 14"	2" x11"	50	50/cs	A6204	

indications:

Pressure ulcers
Leg ulcers
Chronic secreting lesions
I.V. sites
Laparotomy sites

Arthroscopic incisions
Tube sites
Diabetic ulcers
Surgical incisions

contraindications:

Do not use on:
third-degree burns, heavily bleeding wounds, individuals

who are sensitive to or have had an allergic reaction to the dressing or its materials

