

## Cardinal Health™ Negative Pressure Wound Therapy Systems

Seamless care for every patient,  
every step of the way

Simple. Reliable. Effective.

# Cardinal Health™ Negative Pressure Wound Therapy Systems


## Simple. Reliable. Effective.

To help promote healing from hospital to home, Cardinal Health™ has a portfolio of Negative Pressure Wound Therapy (NPWT) devices and disposables to meet your patients' needs across the continuum of care. Our NPWT portfolio provides a variety of options that deliver a market-leading standard of care — without escalating costs.

## Top five reasons to use Cardinal Health™ NPWT

A unique balance between your facility's clinical and financial goals.

- 1 Clinically equivalent technology** to the market leader.  
Our devices have proven the test of time, and have been FDA cleared since 2007.<sup>1</sup>
- 2 Save an average of 30–50%** on your NPWT spend compared to the market leading products.<sup>2</sup>
- 3 We offer easy to use products** that provide a common patient and clinician experience across the continuum of care.
- 4 Our patented Cardinal Health™ Simultaneous Irrigation™ Technology** is the only option that delivers a topical solution into the wound bed at the same time our Negative Pressure Wound Therapy system is applied.<sup>3</sup>
- 5** The Cardinal Health team is **dedicated to assisting you** through the discharge process.





## Device Portfolio

Cardinal Health™ NPWT Products target the needs of patients in different care settings while ensuring consistency across the continuum of care. Our devices share similar designs and user interfaces so the products you choose deliver a seamless experience for you and your patients.

 <p><b>SVED®</b></p>	 <p><b>CATALYST™</b></p>	 <p><b>PRO</b></p>	 <p><b>PRO TO GO</b></p>	 <p><b>ALLY™</b></p>	 <p><b>ALLY TO GO™</b></p>
<p><b>Reusable NPWT systems</b></p>		<p><b>Disposable NPWT systems</b> – TO GO devices are single use –</p>			

CATALYST™, ALLY™ and ALLY TO GO™ include occlusion detection



## Cardinal Health™ SVED® and CATALYST™ Therapy System

For a reliable, reusable NPWT device, choose SVED® or CATALYST™:

Our technology has provided over 1.3 million days of therapy — while offering flexible purchase options to save up to 30% to 50% compared to market-leading devices<sup>2</sup>. These devices offer intuitive design and easy to navigate button controls, reducing complexity for clinicians, patients, and caregivers.

- Small, quiet, light-weight device
- Long battery life (up to 18 hours)
- Visual and audible alarms
- Continuous and intermittent therapy modes
- Can be used with simultaneous wound irrigation
- Available with or without occlusion detection

### Ordering information

Cat. No.	Item
6701132	Cardinal Health™ SVED® NPWT Device
68-1132	Cardinal Health™ CATALYST™ NPWT Device



## Simultaneous Irrigation™ Technology

Some wounds might call for more than just NPWT. With our unique, patented Simultaneous Irrigation™ Technology, you have the option to deliver an irrigant into the wound bed at the same time NPWT is applied.

Our Simultaneous Irrigation™ Technology:

- Reduces bioburden and removes debris<sup>5</sup>
- Moisturizes and cleanses the wound
- Allows continuous NPWT and helps maintain dressing seal integrity
- Can be used with any Cardinal Health™ NPWT device




## Make a best practice even better

Irrigation is already a best practice for complex wound cleansing. Simultaneous Irrigation™ Technology takes it to the next level.




### Ready

Apply NPWT dressing and close irrigation roller clamp completely. Add solution to bag. Cut hole in drape and apply Simultaneous Irrigation™ SpeedConnect™.



### Set

Turn device on and select negative pressure setting and mode. Allow dressing to reach set pressure.



### Flow

Open the clamp and allow irrigation solution to flow. Set the drip rate per physician order.



## Ordering information

Cat. No.	Item
47-6000	Cardinal Health™ NPWT Irrigation Tubing with SpeedConnect™
47-6500	Cardinal Health™ NPWT Irrigation Delivery Set

## Cardinal Health™ PRO and ALLY™ Therapy Systems

Our disposable NPWT systems provide a high-performance device in a lightweight, portable size. The disposable PRO and ALLY™ devices are:

- Designed for up to 42 days with 1,000 hours of use
- Reusable across multiple patients
- Conveniently portable
- Easy to use, with simple buttons
- Equipped with visual and audible alerts
- Available with or without occlusion detection
- Compatible with Simultaneous Irrigation™ Technology
- Created to have up to 24 hours of battery life



### PRO and ALLY™ ordering information

Cat. No.	Item
6708888	Cardinal Health™ PRO NPWT Device
68-8888	Cardinal Health™ ALLY™ NPWT Device

## Cardinal Health™ PRO TO GO and ALLY TO GO™ Kits

Our TO GO Kits are convenient, single-patient use solutions that can help minimize the risk of cross-contamination. Our TO GO Kits:

- Contain everything you need for 10 days of therapy
- Include 10 days of therapy, three dressing kits, two canisters and a carrying case
- Are available with or without occlusion detection
- Can be used as a therapy starter kit or for transition between sites of care



### PRO TO GO and ALLY TO GO™ ordering information

Cat. No.	Item
47-1230	Cardinal Health™ PRO TO GO NPWT Kit (includes one 10 day PRO TO GO NPWT Device, three large dressings and two 300cc canisters)
48-1230	Cardinal Health™ ATG™ NPWT KIT (includes one 10 day ALLY TO GO™ NPWT Device, three large dressings and two 300cc canisters)



## Disposables and Accessories

We offer a comprehensive line of disposable devices and high-quality accessories designed to meet your patients' needs and integrate with your NPWT devices across the continuum of care.



- Black foam dressing kits
- White foam
- SensiSkin™ drape
- Bridging kits
- SpeedConnect™ tubing
- Canisters
- Y connectors
- IV poles adapters
- Carrying cases

### SVED® and PRO disposables and accessories ordering information

Cat. No.	Item
47-1702	Cardinal Health™ NPWT Black Foam Dressing Kit, small
47-1701	Cardinal Health™ NPWT Black Foam Dressing Kit, medium
47-1700	Cardinal Health™ NPWT Black Foam Dressing Kit, large
47-1703	Cardinal Health™ NPWT Black Foam Dressing Kit, XL
47-1702-15	Cardinal Health™ NPWT Black Foam Dressing Kit, small, OD
47-1701-15	Cardinal Health™ NPWT Black Foam Dressing Kit, medium, OD
47-1700-15	Cardinal Health™ NPWT Black Foam Dressing Kit, large, OD
47-4000	Cardinal Health™ NPWT 300 cc Canister with Gel
47-4100	Cardinal Health™ NPWT 300 cc Canister without Gel
47-4500	Cardinal Health™ NPWT 500 cc Canister with Gel
47-2000	Cardinal Health™ NPWT SpeedConnect™ Tubing
47-2500	Cardinal Health™ NPWT Y Connector

### CATALYST™ and ALLY™ disposables and accessories ordering information

Cat. No.	Item
48-1702	Cardinal Health™ NPWT Black Foam Dressing Kit, small, OD
48-1701	Cardinal Health™ NPWT Black Foam Dressing Kit, medium, OD
48-1700	Cardinal Health™ NPWT Black Foam Dressing Kit, large, OD
48-1703	Cardinal Health™ NPWT Black Foam Dressing Kit, XL, OD
48-1702-15	Cardinal Health™ NPWT Black Foam Dressing Kit, small, OD
48-1701-15	Cardinal Health™ NPWT Black Foam Dressing Kit, medium, OD
48-1700-15	Cardinal Health™ NPWT Black Foam Dressing Kit, large, OD
48-4000	Cardinal Health™ NPWT 300 cc Canister with Gel, OD
48-4100	Cardinal Health™ NPWT 300cc Canister without Gel
48-4500	Cardinal Health™ NPWT 500 cc Canister with Gel, OD
48-2000	Cardinal Health™ NPWT SpeedConnect™ Tubing, OD
48-2500	Cardinal Health™ Occlusion Detection Y Connector

### White foam and accessories ordering information

Cat. No.	Item
47-1704	Cardinal Health™ NPWT Bridging Kit
47-1751	Cardinal Health™ NPWT White Foam Dressing, small
47-1755	Cardinal Health™ NPWT White Foam Dressing, large
47-5500	Cardinal Health™ NPWT CATALYST/SVED IV Pole Holder
47-5600	Cardinal Health™ NPWT ALLY/PRO Family IV Pole Holder
47-7000	Cardinal Health™ NPWT Polyurethane Drape
47-7100	Cardinal Health™ NPWT SensiSkin™ Drape
47-9000	Cardinal Health™ NPWT SVED®/CATALYST™ A.C. Power Adapter
47-9100	Cardinal Health™ NPWT ALLY™/PRO Family A.C. Power Adapter
47-9500	Cardinal Health™ NPWT CATALYST™/SVED® Carrying Case
47-9600	Cardinal Health™ NPWT ALLY™/PRO Family Carrying Case



## Supporting the continuum of care

### Discharge providers

Cardinal Health™ NPWT knows how important discharge can be. We offer several simple and convenient options to make the discharge process seamless for you and your patients. Please contact your Cardinal Health NPWT representative to learn more about our discharge providers who can help coordinate continued therapy as your patients transition out of your facility.

### NPWT C.A.R.E. Program

Our Continued Access to Resources and Education (C.A.R.E.) program is designed to make it easy for you to learn more about Cardinal Health™ NPWT and get the on-going support you need. Visit [cardinalhealth.com/NPWTcare](http://cardinalhealth.com/NPWTcare) for the full experience.

### Product support

Reach us 24/7 for product support at 866-484-6798.

### Virtual in-services and training

Virtual staff education is available. We can accommodate virtual training times for the entire nursing team or individual trainings with a Cardinal Health Clinician.

### Instructional videos

Watch and learn about the proper use of our devices, including:

- Device cleaning
- Device-specific operations
- Troubleshooting
- Simultaneous Irrigation™ Technology
- Dressing application
- Bridging and Y Connector Techniques

### Continuing Education courses

Stay up to date and earn CE credit with courses on a range of wound management topics, such as pressure injuries, wound assessment, and advanced NPWT techniques.

Learn more at [cardinalhealth.com/NPWT](http://cardinalhealth.com/NPWT)

### Clinical research summary

To compare clinical equivalence across NPWT devices, Dr. Lawrence Lavery, DPM, MPH at UT Southwestern, performed a randomized clinical study, funded by Cardinal Health. The study found no significant difference in clinical outcomes or adverse events when comparing NPWT devices<sup>1</sup>.

### References

- 1 Davis, K. E., La Fontaine, J., Farrar, D., Oz, O. K., Crisologo, P. A., Berriman, S., & Lavery, L. A. (2020). Randomized clinical study to compare negative pressure wound therapy with simultaneous saline irrigation and traditional negative pressure wound therapy for complex foot infections. *Wound Repair and Regeneration*, 28(1), 97-104. doi:10.1111/wrr.12741. Epub 2019 Jun 27.
- 2 Data on file at Cardinal Health
- 3 Patents: US7837673, US7608066, US7532953, US8142405
- 4 Occlusion detection features dual-lumen tubing which maintains and monitors pressure at the wound site and allows for continuous removal of exudate during NPWT
- 5 Davis KE, Moquin KJ, Lavery LA. The fluid dynamics of simultaneous irrigation with negative pressure wound therapy. *Int Wound J*. 2016;13(4):469-474.

© 2020 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, ESSENTIAL TO CARE, ALLY, CATALYST, KENDALL, SENSISKIN, SIMULTANEOUS IRRIGATION, SPEEDCONNECT and SVED are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other trademarks are the property of their respective owners. Lit No. 2NP20-1263911

