




# Designed for more **placement flexibility** than the PleurX<sup>®</sup> System.

With three placement method options (antegrade, retrograde and over-the-wire), and multiple kit options designed to accommodate customer and patient needs, the Aspira® Drainage System gives you more placement flexibility than the PleurX® Drainage System.




## More consistent pressure during drainage.

The Aspira Drainage System is designed to drain quickly, but without the use of high pressure. Lower, consistent vacuum pressure is preferable because high vacuum pressure may cause patient discomfort.

- Simulated testing has shown that the Aspira Drainage System has a 68% higher average flow rate than the PleurX Drainage System<sup>1,2,4</sup>
- Simulated testing of the Aspira Drainage System demonstrated consistent pressure and flow rate, unlike the PleurX Drainage System, which has high initial vacuum that drops quickly<sup>3</sup>
- The Aspira Drainage System has 67% more fenestrations and a larger valve than the PleurX Drainage System



**ASPIRA Drainage System** 



**72%** larger valve

Simulated Average Flow Rate Comparison		
ASPIRA Drainage System	PleurX Drainage System	
7.46 mL/sec	4.45 mL/sec	

68%
higher
average flow
rate²

PleurX Drainage System

ASPIRA Drainage System

67%
more
fenestrations

- 1. Data on file. Internal diameters measured at the narrowest point. N=10 per device.
- 2. Data on file. Simulated water drainage. N=9 for Aspira Drainage Bag, N=8 for PleurX Drainage Bottle. Bench data may not be represented of clinical outcomes. Different test methods may yield different results.
- 3. Data on file. Simulated water drainage under similar conditions. N=40 for Aspira Drainage Bag. N=8 for PleurX Drainage Bottle. Bench data may not be represented of clinical outcomes. Different test methods may yield different results.
- 4. Data on file.



# Convenient options not offered by the PleurX® System.

#### SIMPLE TO USE

With its simple and easy-to-use design, the Aspira Drainage System is an excellent choice for caregivers and patients.



#### EASY OCCLUSION MANAGEMENT

With the Aspira Drainage System there is no risk of accidentally discharging vacuum. Vacuum is activated with one squeeze of the manual pump.

Replaceable valve allows for catheter repair - rather than catheter exchange.



#### **ASPIRA Drainage System**


7.51 NV Drainage System		
Peritoneal Product Code	Pleural Product Code	Description
4992206	4992506	Catheter Kit 1 Insertion Tray with Catheter
4992207	4992507	Complete Kit 1 Insertion Tray with Catheter, 5 Drainage Kits, 5 Dressing Kits, Voucher for 20 Drainage Kits and 20 Dressing Kits
4992208	4992508	Essential Kit 1 Insertion Tray with Catheter, 20 Drainage Kits, 20 Dressing Kits
4992209	4992509	Starter Kit 1 Insertion Tray with Catheter, 5 Drainage Kits, 5 Dressing Kits

less storage space\*

Twenty Aspira drainage bags and dressings require approximately the same amount of space as a box of just 4 PleurX bottles and dressings reducing storage space and waste.

#### **ASPIRA** Drainage Accessories

Peritoneal Product Code	Description
4991503	Dressing Kit – 5 Kits Included
4992301	1000 mL Drainage Bag – 5 Bags Included
4992305	Luer Adapter
4992306	Repair Kit
4992307	Universal Tubing Adapter



#### Visit myaspira.com for more information.

\* Data on file. PleurX packaging measurements represent five boxes of PleurX Starter Kit containing four 1,000 mL drainage bottles and dressings per box. Aspira packaging measurements represent one box of Aspira Essential Kit containing twenty 1,000 mL drainage bags and dressings per box.

All trademarks are the property of their respective owners.

#### Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



### **ENDOTFK®**

#### Understand. Innovate. Deliver.

Merit Medical Systems, Inc. 1600 West Merit Parkway South Jordan, Utah 84095 1.801.253.1600 1.800.35.MERIT

myaspira.com

Merit Medical Europe, Middle East, & Africa (EMEA)

Amerikalaan 42, 6199 AE Maastricht-Airport The Netherlands

+31 43 358 82 22

Merit Medical Ireland Ltd.

Parkmore Business Park West Galway, Ireland +353 (0) 91 703 733

+61 (0) 3 8373 4943

**Austria** 0800 295 374

Belgium

0800 72 906 (Dutch) 0800 73 172 (Français)

+55 11 4561 - 0062

Canada 1.800.364.4370

+86-10 8561 0788

Denmark 80 88 00 24 Finland 0800 770 586

0800 91 60 30

Germany

0800 182 0871

Hong Kong +852-2207 6888

+91-80-41223376

Ireland (Republic)

1800 553 163

Italy 800 897 005

Luxembourg 8002 25 22

Netherlands 0800 022 81 84

Norway 800 11629

Portugal 308 801 034

+7 495 221 89 02

Singapore +65 67536311

South Korea +8210 2205 0707

(Italiano) +41 22 518 00 35

0800 973 115

+34 911238406

+41 22 518 02 30

(Français) +41 22 518 02 52

Sweden 020 792 445

Switzerland

(Deutsch)