

AIRCAST

■ Open-frame design for ventilation

AirSelect™ Walker Boots disperse pressure evenly with a 3-strap system. Ultra-light, trimmable, heat-moldable, semirigid shell uses an open-frame design for ventilation. Wearers can use integrated inflation system to inflate each aircell individually to maximize comfort and minimize edema. Available in three styles: **Short**, **Standard**, and **Elite**. **Short** offers moderate support for foot injuries. Designed specifically for metatarsal fracture, forefoot and midfoot injury, acute or post-operative use, bunionectomy, and soft tissue injury. **Standard** offers moderate support for lower leg, foot, and ankle injuries. Best for stable fractures of the lower leg, foot, and ankle; severe ankle sprains; and post-operative use. **Elite** is engineered to provide ultimate protection, comfort, and edema control. Designed for stable fractures of the lower leg, foot, and ankle; severe ankle sprains; and post-operative use. Enclosed toe cover provides additional protection. Fits left or right.



SIZE	MENS SHOE SIZE	WMNS SHOE SIZE	SHORT	STANDARD	ELITE	QTY
X-SMALL	UP TO 4	UP TO 5	#2970006836	#2970006841	#2970006846	EACH
SMALL	4-7	5-8	#2970006837	#2970006842	#2970006847	EACH
MEDIUM	7-10	8-11	#2970006838	#2970006843	#2970006848	EACH
LARGE	10-13	11-15	#2970006839	#2970006844	#2970006849	EACH
X-LARGE	13+	15+	#2970006840	#2970006845	#2970006850	EACH

IMC • California residents, see left page.

Accessories

Two-layer oversole adjusts height



■ Compensates for temporary leg-length discrepancy

EVENUP™ Orthotic Shoe Lifts add ½" to 1" lift to compensate for a temporary leg-length discrepancy such as when wearing a walking cast or fracture boot. Easily attaches to shoe of nonaffected limb for improved balance and gait, and reduced body strain. Rubber and hook-and-loop straps hold it in place. Includes adjustable two-layer nonskid rubber outsole. **Pro Pack** includes five of each size. Fits left or right.

Sizing: Small: Sole measurement between 10" and 11¼"; Wmns 6-8½, Mens 7-8½. Medium: Sole measurement between 11¼" and 12¼"; Wmns and Mens 9-11½. Large: Sole measurement between 13" and 14½"; Wmns and Mens 12-14.*

#2970005920	SMALL	EACH
#2970005921	MEDIUM	EACH
#2970005922	LARGE	EACH
#2970005923	PRO PACK	15/CASE

IMC • California residents, see left page.

*It is generally best to choose an EVENUp a little too small rather than a little too big. This means if your shoe measurement is between two sizes, like 11¼"-11½", or measures 11½", it is best to go with the small. Also note that shoe size best correlates with size, as well.

■ Peel-away progressive adjustment

Achilles Boot Wedges fit into most closed-heel walker boots, turning any walker boot into an Achilles Tendon Walker Boot. Hindfoot wedge has peel-away layers that allow progressive adjustment of plantarflexion. Foot is initially positioned in 30° of plantarflexion. As recovery progresses, decrease plantarflexion in ½° increments by peeling a layer away from the base of wedge. Fits left or right.



Sizing: Medium, Wmns up to 11, Mens up to 10; Large, Wmns 11+, Mens 10+.

#2970006242	MEDIUM	EACH
#2970006243	LARGE	EACH

IMC • California residents, see left page.

DARCO


■ Fits most standard walker boots

PegAssist™ Walker Off-Loading Insole System Multi-indication removable peg chassis effectively off-loads the plantar aspect of the foot after surgery or when wounds are present. Reduces pressure by as much as 60% while still allowing wearer to remain ambulatory. Stabilizer board has a low-profile molded hook strip and adheres to liner to prevent shifting while walking. Universal system allows the system to fit most standard walker boots.



Removable pegs for a custom fit

Sizing: X-Small, Wmns 3½-6½, Mens 2-4½; Small, Wmns 6½-9, Mens 5-7½; Medium, Wmns 9½-12, Mens 8-10½; Large, Wmns 12½-15, Mens 11-13; X-Large, Wmns 15+, Mens 14+.



#2970007251	X-SMALL	EACH
#2970007250	SMALL	EACH
#2970007252	MEDIUM	EACH
#2970007249	LARGE	EACH
#2970007253	X-LARGE	EACH

IMC • California residents, see left page.