


Wheelchair Accessories



Keep fingers out of the spokes



Wheelchair Spoke Covers

are essential for any patient with impaired sensation or alertness. Fabric-covered foam disks cover spokes so fingers cannot get caught between them. Attach in seconds with internal hook-and-loop fasteners.

May ship in Black or Blue

Sizing: Measure diameter of hand rim. Small, 20" diam.; Large, 22" diam.

#2970009888 SMALL PAIR

#2970009887 LARGE PAIR

MD • **California residents, see below.**



Eliminate "pinch points"



Use the palms of hands to push wheel

Shove-a-Lugs® Wheelchair Grip Attachments

Grips work by enabling user to push the wheel with the palms of the hands, rather than gripping the thin rim at an awkward angle with a pinch grip. Important for people with arthritic, painful, or weakened hands. No gripping is required—an ergonomic solution to a non-ergonomic mechanism.

Contains: 12 Shove-a-Lug grip attachments, 24 hex screws, 1 hex wrench, 12 adapter rings

#2970008933 WHEELCHAIR GRIP ATTACHMENT PACK

California residents, see below.



Improves grip

Wheel-Ease Rim Covers This alternative to gloves protects the hands from calluses and provides a comfortable handhold. Foam padding and contoured grooves improve hand grip and aid in propulsion. Helps protect against injury by shielding hands from hard metal rims. Easy, tool-free installation. Black.

#2970009894 WHEEL-EASE RIM COVER PAIR

California residents, see below.



Maximum impact protection

Fingerless Gloves and Inserts LYCRA® base with viscoelastic polymer pad over palm and between thumb and index finger. Webbing is covered with grain leather. Maximum impact protection. Optional Inserts sold separately.

Black

Sizing: Measure hand circumference across knuckles at the base of the fingers. Small, 7"-8"; Medium, 8"-9"; or Large, 9"-10".

#2970009821 GLOVES, SMALL PAIR

#2970009820 GLOVES, MEDIUM PAIR

#2970009819 GLOVES, LARGE PAIR

#2970009818 INSERTS, SMALL PAIR

#2970009817 INSERTS, MEDIUM PAIR

#2970009816 INSERTS, LARGE PAIR

California residents, see below.



Durable foam-padded leather

Heavy-Duty Wheelchair Gloves protect palm while leaving fingers free. Designed to be easy to put on and remove. High-quality materials provide comfort and performance.

Foam-padded Leather on wrist and palm • Breathable spandex on back of hand • Black

Sizing: Measure circumference of MCP joints. Small/Medium, 7"-8"; Large/X-Large, 8 1/2"-10"

#2970010955 SMALL/MEDIUM PAIR

#2970010988 LARGE/X-LARGE PAIR

California residents, see below.

Smoker's Aprons



Heavy-Duty Smoker's Aprons

Fire-resistant apron covers wheelchair user's lap and chest to help prevent smoking materials from falling between them and wheelchair. Ties around neck, wheelchair arms, and just below seat. Can be folded safely for easy storage. **Economy Smoker's Aprons** are a lightweight, more economical version and cannot not be folded due to fiberglass content.

Flame-retardant fabric meets NFPA 701 • **Heavy-Duty:** Ships folded • Stores folded • Heavy-duty material • No fiberglass content • Grey • **Economy:** Store only by hanging due to fiberglass content • Flame-retardant fabric meets NFPA 701 • White

#2970011515 HEAVY-DUTY EACH

#2970011514 ECONOMY EACH

California residents, see below.