

Wrist Immobilizers

Adjustable dorsal and palmar stays

FREEDOM® comfort™ Wrist Supports offer optimum support and immobilization for wrist joint. Adjustable dorsal and palmar metal stays for stabilization of wrist joint. Narrow web space helps prevent pinching. Contoured along distal palmar crease to allow 90° MP flexion. Open-cell foam padding with breathable polyester laminate against skin. Hook-and-loop D-ring closures. Not made with natural rubber latex.

Open-cell foam pad • Breathable polyester laminate liner

Sizing: Measure circumference of hand at MP joints.
See chart below.

SIZE	HAND CIRC.	LEFT	RIGHT	QTY
SMALL	7¼"-8¼"	#2970010195	#2970010191	EACH
MEDIUM	8¼"-9¼"	#2970010196	#2970010192	EACH
LARGE	9¼"-10¼"	#2970010197	#2970010193	EACH
X-LARGE	10¼"-12"	#2970010198	#2970010194	EACH

MD • **CA** California residents, see left page.



Centered palmar stay for sturdy, rigid support



After-cast immobilization

Wrist Immobilizers Constructed of grey perforated vinyl laminated to a soft flannel lining with malleable aluminum palmar and dorsal stays **Short (8")** Use after cast removal for wrist fractures or for mild-to-severe wrist sprains/strains and carpal tunnel syndrome. **Long (10")** Use after cast removal for distal forearm fractures. Additional wrist strap for greater support.

Sizing: Measure wrist circumference. See charts below.

SIZE	HAND CIRC.	LEFT	RIGHT	QTY
SHORT, 8"				
SMALL	5½"-6½"	#2970003470	#2970003469	EACH
MEDIUM	6½"-7¼"	#2970003472	#2970003471	EACH
LARGE	7¼"-8"	#2970003474	#2970003473	EACH
X-LARGE	8"-8¾"	—	—	EACH

CA California residents, see left page.

SIZE	HAND CIRC.	LEFT	RIGHT	QTY
LONG, 10"				
SMALL	5½"-6½"	—	—	EACH
MEDIUM	6½"-7¼"	#2970003094	#2970003085	EACH
LARGE	7¼"-8"	#2970003128	#2970003111	EACH
X-LARGE	8"-8¾"	#2970003182	#2970003147	EACH

CA California residents, see left page.



Short, 8"




Long, 10" has extra wrist strap for greater support

Use the FREEDOM® Long for:

- ✓ All-day comfort
- ✓ Relief from carpal tunnel syndrome and tenosynovitis
- ✓ Protecting the wrist after cast removal
- ✓ Alleviating symptoms of arthritis affecting the wrist

Long Elastic Wrist Supports The longer design is perfect for all-day support. Long forearm length optimally distributes pressure and hand weight without restricting hand function. Offers an additional benefit to those with weak or arthritic wrists.



Long Elastic

SIZING (MEASURE WRIST CIRC.)		USA SHORT (RED/BLUE)		SHORT ELASTIC (BEIGE)		LONG ELASTIC (BEIGE)		QTY
		LEFT	RIGHT	LEFT	RIGHT	LEFT	RIGHT	
XX-SMALL	3½"-4½"	—	—	—	—	#2970003593	#2970003592	EACH
X-SMALL	4½"-5½"	#2970002261	#2970002262	#2970003607	#2970003606	#2970003595	#2970003594	EACH
SMALL	5½"-6½"	#2970002267	#2970002268	#2970003609	#2970003608	#2970003597	#2970003596	EACH
MEDIUM	6½"-7¼"	#2970002264	#2970002266	#2970003611	#2970003610	#2970003599	#2970003598	EACH
LARGE	7¼"-8"	#2970002263	#2970002265	#2970003613	#2970003612	#2970003601	#2970003600	EACH
X-LARGE	8"-8¾"	—	—	—	—	#2970003603	#2970003602	EACH

MD • **CA** California residents, see left page.

LOG IN FOR PRICING at AliMed.com