

## Standard Percutaneous Endoscopic Gastrostomy (PEG) Kit

## **Product Features**

# High-Grade Silicone Dome-Bolstered PEG Tube

- Designed for performance and consistent flow
- > Facilitates traction removal

## Universal 1.9 mm O.D. x 240 cm Retrieval Snare

- Traverses scopes with a 2.0 mm or larger working channel
- ➤ Universal-sized snare, intended to eliminate the cost of a specialty snare

## **Double-Tier Tray Packaging**

 Kit organization promotes procedural efficiency

## **External Bolsters**

➤ Two external bolster options: accommodates physician and patient needs: Round; T-Bar

## **Multiple Sizes and Accessories**

- Multiple French sizes for broad patient coverage
- > Separate feeding ports
- C-Clamp designed to help reduce risk of tube reflux and leaking

#### Standard PEG Kit

Non-ENFit Codes	ENFit Codes	Description	Units
M005 <b>6820</b> 1	M005 <b>0904</b> 1	20Fr Standard PEG Kit, Pull with Lidocaine Ampule	Box 2
M005 <b>6821</b> 1	M005 <b>0905</b> 1	20Fr Standard PEG Kit, Push with Lidocaine Ampule	Box 2
M005 <b>6824</b> 1	M005 <b>0906</b> 1	24Fr Standard PEG Kit, Pull with Lidocaine Ampule	Box 2
M005 <b>6825</b> 1	M005 <b>0907</b> 1	24Fr Standard PEG Kit, Push with Lidocaine Ampule	Box 2

## **PEG Accessories**

Non-ENFit Codes	ENFit Codes	Description	Units
M005 <b>8064</b> 1	M005 <b>0126</b> 1	20Fr Y-Ports	Box 10
M005 <b>8065</b> 1	M005 <b>0127</b> 1	24Fr Y-Ports	Box 10
M005 <b>8073</b> 1	N/A	20Fr Round Bolsters	Box 10
M005 <b>8078</b> 1	N/A	C-Clamps (Fits 20Fr and 24Fr)	Box 10

## **EndoVive**

## **Enteral Access Devices**

## Safety Percutaneous Endoscopic Gastrostomy (PEG) Kit

#### **Product Features**

## **Deluxe Safety PEG Procedure Kit**

- ➤ High-grade silicone dome-bolstered PEG tube
- ➤ ChloraPrep® Skin Preparation
- ➤ Universal 1.9 mm O.D. x 240 cm retrieval snare
- Double-tier tray packaging
- ➤ Two external bolsters
- ➤ Multiple sizes and accessories
- ➤ Available for push and pull methods

## **Safety Trocar Catheter**

- > Needle tip protector designed to engage when the needle stylet is removed from the needle cannula
- Designed to allow re-activation by replacing the needle stylet back into the needle cannula

#### Safety-Shielded Scalpel


➤ Thumb-activated safety mechanism retracts spring-loaded No. 11 surgical blade into handle

## SafetyGlide® Shielded **Injection Needle**

- > 22 gauge x 1.5" needle designed for injection into abdominal wall
- > Shielding mechanism designed for one-finger activation

#### **PEG Accessories**

- > Y-Ports
- ➤ Round Bolsters
- ➤ T-Bar Bolster
- ➤ C-Clamps



#### Standard PEG Kit

Non-ENFit Codes	ENFit Codes	Description	Units
M005 <b>6820</b> 1	M005 <b>0904</b> 1	20Fr Standard PEG Kit, Pull with Lidocaine Ampule	Box 2
M005 <b>6821</b> 1	M005 <b>0905</b> 1	20Fr Standard PEG Kit, Push with Lidocaine Ampule	Box 2
M005 <b>6824</b> 1	M005 <b>0906</b> 1	24Fr Standard PEG Kit, Pull with Lidocaine Ampule	Box 2
M005 <b>6825</b> 1	M005 <b>0907</b> 1	24Fr Standard PEG Kit, Push with Lidocaine Ampule	Box 2

#### **PEG** Accessories

Non-ENFit Codes	ENFit Codes	Description	Units
M005 <b>8064</b> 1	M005 <b>0126</b> 1	20Fr Y-Ports	Box 10
M005 <b>8065</b> 1	M005 <b>0127</b> 1	24Fr Y-Ports	Box 10
M005 <b>8073</b> 1	N/A	20Fr Round Bolsters	Box 10
M005 <b>8078</b> 1	N/A	C-Clamps (Fits 20Fr and 24Fr)	Box 10