Self-Cath® Plus Coudé Tip

This coudé tip intermittent catheter has fire-polished eyelets designed for comfort and a siliconized surface created for smooth insertion. The catheter has a hydrophilic coating that activates immediately upon exposure to water. The blue Guide Stripe® along the top of the catheter allows monitoring of the coudé tip throughout the catheterization process. This catheter is not made with natural rubber latex.



Self-Cath® Plus Tapered Coudé Tip with Guide Stripe®

Code	Size	Units	Description	HCPCS*
4608	• 8 FR	30	16" length, funnel end	A4352
4610	● 10 FR	30	16" length, funnel end	A4352
4612	O 12 FR	30	16" length, funnel end	A4352
4614	• 14 FR	30	16" length, funnel end	A4352



Self-Cath® Plus Olive Coudé Tip with Guide Stripe®

Code	Size	Units	Description	HCPCS*
4808	• 8 FR	30	16" length, funnel end	A4352
4810	• 10 FR	30	16" length, funnel end	A4352
4812	O 12 FR	30	16" length, funnel end	A4352
4814	• 14 FR	30	16" length, funnel end	A4352
4816	16 FR	30	16" length, funnel end	A4352

Self-Cath® Coudé Tip

Self-Cath coudé tip intermittent catheters have fire-polished eyelets designed for comfort, a siliconized surface created for smooth insertion. The blue Guide Stripe® along the top of the catheter allows monitoring of the coudé tip throughout the catheterization process. They are available in a variety of options. This catheter is not made with natural rubber latex.

Self-Cath® Olive Coudé Tip with Guide Stripe®



Self-Cath® Tapered Coudé Tip with Guide Stripe®

Code	Size	Units	Description	HCPCS*
608	• 8 FR	30	16" length, funnel end	A4352
610	● 10 FR	30	16" length, funnel end	A4352
612	O 12 FR	30	16" length, funnel end	A4352
614	• 14 FR	30	16" length, funnel end	A4352

WARNING: This product can expose you to Di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov