



Medication preparation

Medication delivery



Sharps disposal

Cardinal Health™ Monoject™ Magellan™ Hypodermic Safety Needles

Clinically preferred

Safety needles designed to protect



Your facility depends on solutions that offer greater protection against needlesticks and bloodborne pathogens

Because more than 40% of needlestick injuries happen after use,¹ we designed Cardinal Health™ Monoject™ Magellan™ Safety Needles to keep patients and clinicians safe at every step. Magellan™ Safety Needles are engineered to comply with OSHA regulations and help you standardize solutions across every job with one safety technology.

Choose Monoject™ Magellan™ Hypodermic Safety Needles

Magellan™ Products are the clinically preferred premium option for safety needles and syringes. Intuitive, easy-to-use designs offer protection with just the push of a finger, improving clinician comfort and confidence with less follow-up training.


- 1 Robust safety shield provides coverage of entire cannula and is not easily removed
- 2 One-handed activation, three ways
- 3 Keeps hand behind needle
- 4 Audible activation cue
- 5 Safety shield stays in place, even if the needle bends or breaks
- 6 Design limits blood or fluid splatter

Every Monoject™ Magellan™ Hypodermic Safety Needle does not require a change in protocol and may be activated using the thumb, finger or flat surface.

One-handed activation



Thumb



Finger



Flat surface

Cardinal Health™ Monoject™ Magellan™ Hypodermic Safety Needles

Ordering information

The safety needle and syringe combinations come preattached and safety-engineered to protect clinicians from accidental needlestick injuries.

Features:

- Sterile, single use and individually packaged
- Not manufactured with DEHP or natural rubber latex



Monoject™ Magellan™ Hypodermic Safety Needle



Monoject™ Magellan™ 1 mL Syringe with Hypodermic Safety Needle

Monoject™ Magellan™ Hypodermic Safety Needles

Order code	Description	Needle gauge	Needle length	Box qty.	Case qty.
8881850810	18 G x 1 in. (1.2 x 25 mm)	18 G	1 in.	50	500
8881850815	18 G x 1½ in. (1.2 x 38 mm)	18 G	1½ in.	50	500
8881850910	19 G x 1 in. (1.1 x 25 mm)	19 G	1 in.	50	500
8881850915	19 G x 1½ in. (1.1 x 38 mm)	19 G	1½ in.	50	500
8881850010	20 G x 1 in. (0.9 x 25 mm)	20 G	1 in.	50	500
8881850015	20 G x 1½ in. (0.9 x 38 mm)	20 G	1½ in.	50	500
8881850110	21 G x 1 in. (0.8 x 25 mm)	21 G	1 in.	50	500
8881850115	21 G x 1½ in. (0.8 x 38 mm)	21 G	1½ in.	50	500
8881850210	22 G x 1 in. (0.7 x 25 mm)	22 G	1 in.	50	500
8881850215	22 G x 1½ in. (0.7 x 38 mm)	22 G	1½ in.	50	500
8881850358	23 G x ¾ in. (0.6 x 16 mm)	23 G	¾ in.	50	500
8881850310	23 G x 1 in. (0.6 x 25 mm)	23 G	1 in.	50	500
8881850558	25 G x ¾ in. (0.5 x 16 mm)	25 G	¾ in.	50	500
8881850510	25 G x 1 in. (0.5 x 25 mm)	25 G	1 in.	50	500
8881850515	25 G x 1½ in. (0.5 x 38 mm)	25 G	1½ in.	50	500

Monoject™ Magellan™ Hypodermic Safety Needles with Syringes

Order code	Description	Syringe volume	Needle gauge	Needle length	Box qty.	Case qty.
8881811558	1 mL 25 G x ¾ in. (0.5 x 16 mm)	1 mL	25 G	¾ in.	50	200
8881811510	1 mL 25 G x 1 in. (0.5 x 25 mm)	1 mL	25 G	1 in.	50	200
8881811310	1 mL 23 G x 1 in. (0.6 x 25 mm)	1 mL	23 G	1 in.	50	200
8881833810	3 mL 18 G x 1 in. (1.2 x 25 mm)	3 mL	18 G	1 in.	50	400
8881833010	3 mL 20 G x 1 in. (0.9 x 25 mm)	3 mL	20 G	1 in.	50	400
8881833015	3 mL 20 G x 1½ in. (0.9 x 38 mm)	3 mL	20 G	1½ in.	50	400
8881833110	3 mL 21 G x 1 in. (0.8 x 25 mm)	3 mL	21 G	1 in.	50	400
8881833115	3 mL 21 G x 1½ in. (0.8 x 38 mm)	3 mL	21 G	1½ in.	50	400
8881833210	3 mL 22 G x 1 in. (0.7 x 25 mm)	3 mL	22 G	1 in.	50	400
8881833215	3 mL 22 G x 1½ in. (0.7 x 38 mm)	3 mL	22 G	1½ in.	50	400
8881833310	3 mL 23 G x 1 in. (0.6 x 25 mm)	3 mL	23 G	1 in.	50	400
8881833510	3 mL 25 G x 1 in. (0.5 x 25 mm)	3 mL	25 G	1 in.	50	400
8881833558	3 mL 25 G x ¾ in. (0.5 x 16 mm)	3 mL	25 G	¾ in.	50	400

For more information or to order Monoject™ Sharps Safety Products, contact your sales representative, call **800.964.5227** or visit **cardinalhealth.com/monoject**

All hypodermic safety needles and syringes are ISO color coded.

18 G Pink	19 G Cream	20 G Yellow	21 G Green	22 G Black
23 G Blue	25 G Orange	27 G Gray	28 G Teal	30 G Dark Yellow

Access the Monoject™ Magellan™ Hypodermic Safety Needle instructions for use at **cardinalhealth.com/monoject/IFU**

References: 1. NIOSH. Alert: preventing needlestick injuries in health care settings. Cincinnati, OH: US Dept of Health and Human Services; 1999. <https://www.cdc.gov/niosh/docs/2000-108/pdfs/2000-108.pdf?id=10.26616/NIOSH-PUB2000108>.

© 2021 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, MONOJECT and MAGELLAN are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other trademarks are the property of their respective owners. Patent cardinalhealth.com/patents. Lit. No. 2GM21-1495637-03 (07/2021)

