

96002 Spec Sheet

CREATED BY

Guest

This publication was constructed for you using the Brochure Builder application
available at the Smiths Medical website. Visit us to create your own.

www.smiths-medical.com

Saf-T Holder® Devices

JELCO®

Saf-T Holder® devices come with the back-end needle and multi-sample Luer permanently attached to the holder to minimize potential exposures and holder re-use. Available separately or with Saf-T Wing® Blood Collection & Infusion Sets, Saf-T Holder® devices also provide easy needless transfer after a syringe draw.



ORDERING INFORMATION

Box quantity: 50
Minimum order quantity: 200
Shelf life: 5 years

Please see the *Instructions for Use* for a complete listing of the indications, contraindications, warnings and precautions

PRODUCT BENEFITS

- ▮ Back-end needle and Luer permanently attached to the holder to assure compliance with OSHA directives.
- ▮ Provide safe needleless transfer after a syringe draw.
- ▮ Available separately or with Saf-T Wing® Blood Collection Sets.

PRODUCT SPECIFICATIONS

Saf-T Holder® Devices

Reorder Code Description

96000	Saf-T Holder® with male Luer adapter. A case contains 4 boxes of 50.
96000s	Short Saf-T Holder® with male Luer adapter. A case contains 4 boxes of 50.
96002	Saf-T Holder® with female Luer adapter. A case contains 4 boxes of 50.
96002s	Short Saf-T Holder® with female Luer adapter. A case contains 4 boxes of 50.
96004	Saf-T Holder® Blood Culture device with male Luer adapter. A case contains 1 box of 50.
96005	Saf-T Holder® Blood Culture device with female Luer adapter. A case contains 1 box of 50.
96007	Pediatric tube insert. A case contains 10 bags of 10. (Non-Sterile)
96009	Single Use Holder 2 Bags of 500. (check for regional availability)