


PHMB ANTIMICROBIAL FOAM DRESSINGS

Contains PHMB*, a substance which kills and inhibits the growth of bacteria and yeast within the dressing, and is effective against a broad spectrum of microorganisms (gram+, gram- and yeast) that are frequently associated with the bacterial colonization and infection of wounds.

*Polyhexamethylene Biguanide

GEMCORE
360°



Border



Non Border

KEY FEATURES/BENEFITS:

✓ RAPID & POWERFUL ANTIMICROBIAL RELEASE

Total eradication of MRSE, MRSA, VRE, *S.pyogenes*, *S.epidermidis*, *E.coli*, *C.albicans*, *S.aureus*, *Paeruginosa* within 6 hours.

✓ CONFORMABLE HYDROPHILIC FOAM PAD

Impregnated with PHMB, the foam provides sustained antimicrobial action and helps prevent recolonization of the wound bed.

✓ EXCELLENT FLUID HANDLING CAPABILITIES

Excellent absorption rate (31 g/100 cm²/24 hr) helps reduce the number of dressing changes, provides increase wear time (up to 7 days) and reduces the risk of maceration.

✓ PERFORATED WOUND CONTACT LAYER

Clear perforated film prevents the dressing from adhering to the wound bed.

✓ BEVELED EDGES

Ensures the dressing stays in place to prevent accidental removal by shear between clothes and bedding.


✓ WATERPROOF AND BACTERIAL BARRIER

Protects both the patients and healthcare professional from cross contamination. Keeps bacteria and other contaminants away from the wound. Allows the patient to wash or shower.

INDICATIONS:

PHMB foam dressings are indicated for moderate to heavily exuding chronic and acute wounds that are critically colonized or are at risk of infection. The dressings may be used during the healing process on:

- Post-surgical incisions
- Pressure sores
- Venous stasis ulcers
- Diabetic ulcers
- Donor sites
- Abrasions
- Lacerations
- 1st and 2nd degree burns
- Dermatologic disorders
- Other wounds inflicted by trauma
- As a secondary dressing or cover dressing for packed wounds



MODERATE TO HEAVY
EXUDATING WOUNDS



WOUND
TYPE

ORDERING INFORMATION:

Item #	Size	UOM
@ ACRYLIC ADHESIVE BORDER		
AM-WC-PHMBFB33	3" x 3" (7.5 cm x 7.5 cm) Pad	10/box
AM-WC-PHMBFB44	4" x 4" (10 cm x 10 cm) Pad	10/box
□ NON-BORDER, NON-ADHESIVE		
AM-WC-PHMBFNB44	4" x 4" (10 cm x 10 cm) Pad	10/box
AM-WC-PHMBFNB66	6" x 6" (15 cm x 15 cm) Pad	10/box