

## Alternating Pressure Systems

# Cirrus Alternating Pressure Pump and Pad System

The Cirrus Series Pressure Pump and Pad System is designed with advanced technology that alternately inflates and deflates air cells in the pressure pad. The small compact pump conveniently attaches to the bed and is available with fixed or adjustable pressure. The Cirrus Series features an output of four liters per minute with a five minute cycle time and a one year pump warranty.



APP-CFXI



APP-CAL1 Pictured with Retail Box

## Product Features

- Small, compact design with composite hooks located on back of pump
- Pump produces 4 LPM airflow with five minute cycle time
- Bubble cell design delivers comfort while providing superior therapeutic treatment
- Pad is made of long lasting durable heavy gauge premium vinyl
- Variable pressure setting for optimum comfort and support (APP-CAL I model only)
- Quiet pump alternately inflates and deflates the air cells
- 90 day warranty on mattress pad
- One year warranty on pump

## Product Specifications

- Power: 120 VAC
- Input Frequency: 60 Hz
- Pressure Range: 30-60 mmHg (APP-CAL I)
- Air Volume: 4 liters per minute
- Cycle Time: 5 minutes
- Pump Weight: 2.4 lb.
- Pad Weight: 4 lb. (APP-BPD), 5 lb. (APP-DLXP)
- Pump Dimensions: 6.5" x 6.1" x 3.2"
- Pad Dimensions: 78" x 34" with tubing
- Cell Depth: 2.5"
- Weight Capacity: 275 lb. on bubble pad

| Item #    | Description                                        | Pad Type              | Qty  |
|-----------|----------------------------------------------------|-----------------------|------|
| APP-CFXI  | Cirrus I Fixed Pressure Pump and Pad System        | Basic                 | each |
| APP-CAL I | Cirrus II Alternating Pressure Pump and Pad System | Deluxe with end flaps | each |

Whenever bed rails are used as part of any hospital or home care bed system, the potential for bed rail entrapment exists. To reduce the risk, a patient must always be assessed regarding their risk of entrapment before any product is used. Routine patient monitoring is also required once a product is in use. To learn more about the risk of bed rail entrapment, contact your Roscoe Medical sales representative for a copy of our "Bed Rail Safety Guide" or visit US FDA's website at <http://www.fda.gov> and search for "Bed Rail Entrapment."

## Cirrus Series Replacement Parts

| Item #   | Description                                | Model                   | Qty  |
|----------|--------------------------------------------|-------------------------|------|
| APP-BPD  | APP Replacement Pad, Basic                 | APP-CFX I and APP-CAL I | each |
| APP-DLXP | APP Replacement Pad, Deluxe with End Flaps | APP-CAL I and APP-CAL I | each |
| 90037    | Bubble Repair Kit                          | APP-CFX I and APP-CAL I | each |
| 90047    | APP Pump Fitting                           | APP-CFX I and APP-CAL I | each |
| 90043    | Clear Tubing                               | APP-CFX I and APP-CAL I | each |
| 90049    | Power Cord, 10-Foot                        | APP-CFX I and APP-CAL I | each |
| 90051    | Mounting Bracket                           | APP-CFX I and APP-CAL I | each |
| 90055    | Power Switch                               | APP-CFX I and APP-CAL I | each |



**To Order:**  
(800) 3-ROSCOE  
fax: (440) 572-4261  
roscoemedical.com