## extriCARE® NPWT Dressings

## Distributor Sell Sheet





# **Bandage Dressing**

### Features and Benefits:

- Anatomically fitted bandage dressings available in various sizes for different wound types
- Non-adhesive contact layer ensures a tight seal and allows negative pressure to reach the wound bed for optimal healing
- Patented flat, perforated drainage tube ensures even distribution of negative pressure
- Connects via tubing to the ExtriCARE canister for effective removal of exudate, maintaining a cleaner, healthier wound site

HCPCS Code: A6550

# Foam Dressing

### Features and Benefits:

- 2 step drape and easy to follow instructions
- Can customize tubing length to tailor to patient needs
- Low profile soft dome port pad
- High tensile strength foam to reduce flaking when cut to any wound size
- Large, one size fits all kit, (small and medium sizes are available if needed)
- Included ruler sticker with each kit to assist with wound Measurement



EC-Foam-L-E

## Enhancing Wound Care with ExtriCARE®

The ExtriCARE® Hybrid hNPWT technology, allows the use of both traditional foam dressing kits and anatomically fitted bandage dressings with any of the ExtriCARE® NPWT pumps included the ExtriCARE® 1000, 3000, and 3600. This versatility promotes effective wound healing by providing tailored solutions for various wound types.



### **Bandage Dressing Specifications**

Mini: Pad 6 cm x 5 cm (2.36 in x 1.97 in) Small: Pad 10.2 cm x 7.4 cm (4.02 in x 2.91 in) Large: Pad 15.2 cm x 7 cm (5.98 in x 2.76 in) X-Large: Pad 22 cm x 6.5 cm (8.66 in x 2.60 in)

**Tubing Length:** 1.5m (4.92ft)

#### **Dressing Kit Includes:**

- Anatomically fitted bandage dressing
- Connection tubing to ExtriCARE® pump
- · Ruler on packaging Instructions for use

## Foam Dressing Specifications

Size: Large (one size fits all), Medium, Small

Small: 10cm x 7.5cm (3.94in x 2.95in) Medium: 20cm x 13cm (7.87in x 5.11in) Large: 25cm x 16cm (9.84in x 6.30in)

Tubing Length: 2 m (6.56 ft)

#### **Dressing Kit Includes:**

- Anatomically fitted bandage dressing
  Connection tubing to ExtriCARE® pump
- Ruler on packaging
- Instructions for use

· Adhesive strips



Alleva Medical Ltd

Suite M-O. 12th Floor, Kings Wing Plaza 2, 1 On Kwan Street, Shatin, Hong Kong Tel: +852 3166 7200 www.allevamedical.com

Distributed in the USA exclusively by

ExtricareUSA LLC

Suite 205, 2460 General Armistead Avenue, Norristown, PA 19403 Tel: 877-312-NPWT info@extricareusa.com

All products are manufactured under US FDA, ISO 13485, and Japan PMDA registered facility of Alleva Medical Ltd.

Alleva Medical®, extriCARE®, Vide®, Maxson® are Trademarks and/or Registered Trademarks of Alleva Medical Ltd., its subsidiaries, and/or affiliates in the United States, and China and other countries. Copyright 2023. All dimensions, weights, colours and contents may vary.



