

Primer® Modified Unna Boot

Primer Modified Unna Boot is a non-raveling thread locked gauze bandage that has been evenly impregnated with a non-hardening paste of zinc oxide, available with or without calamine. The **Primer Modified Unna Boot** is for use in a compression dressing for the lower leg. Ideal for managing extremities such as varicose ulcers, requiring protection and support.



CODE	SIZE	AMT/BOX	AMT/CASE	HCPCS
GL3001	3" x 10 yds (7.6cm x 9m)	N/A	12	A6456
GL4001	4" x 10 yds (10cm x 9m)	N/A	12	A6456
GL3001C	3" x 10 yds w/Calamine	N/A	12	A6456
GL4001C	4" x 10 yds w/Calamine	N/A	12	A6456

Unna Pak: Primer® Modified Unna Boot and Duban Self-Adherent Bandage

The **Unna Pak** is for the ambulatory management of venous stasis leg ulcers. A convenient system containing two bandages that together form a two-layer Unna Boot. The **Una-Pak** has all the advantages of the **Primer Modified Unna Boot** with the additional benefit of a flesh colored secondary dressing that is lightweight and easy to apply.



CODE	SIZE	AMT/BOX	AMT/CASE	HCPCS
GL2003	3" (7.6cm) Paste Bandage	N/A	12	A4649
GL2004	4" (10cm) Paste Bandage	N/A	12	A4649