

# New Image Two-Piece Pouching System

## Convex Skin Barriers with Integrated Floating Flange



**1140x** (CUT-TO-FIT)

- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border

**Box of 5**



**HCPCS: A4407**

Stock No	Flange Size	For Stomas Up To
11402	1-3/4" (44 mm)	1" (25 mm)
11403	2-1/4" (57 mm)	1-1/2" (38 mm)
11404	2-3/4" (70 mm)	2" (51 mm)



**1530x** (CUT-TO-FIT)

- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ Without Tape Border

**Box of 5**



**HCPCS for 15302, 15303: A4407**

**HCPCS for 15304: A4408**

Stock No	Flange Size	For Stomas Up To
15302	1-3/4" (44 mm)	1" (25 mm)
15303	2-1/4" (57 mm)	1-1/2" (38 mm)
15304	2-3/4" (70 mm)	2" (51 mm)



**1480x** (CUT-TO-FIT)

- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border

**Box of 5**



**HCPCS: A4407**

Stock No	Flange Size	For Stomas Up To
14802	1-3/4" (44 mm)	1" (25 mm)
14803	2-1/4" (57 mm)	1-1/2" (38 mm)
14804	2-3/4" (70 mm)	2" (51 mm)



**1580x** (CUT-TO-FIT)

- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ Without Tape Border

**Box of 5**



**HCPCS for 15802 & 15803: A4407**


**HCPCS for 15804: A4408**

Stock No	Flange Size	For Stomas Up To
15802	1-3/4" (44 mm)	1" (25 mm)
15803	2-1/4" (57 mm)	1-1/2" (38 mm)
15804	2-3/4" (70 mm)	2" (51 mm)



**DID YOU KNOW...**

The ***CeraPlus*** skin barrier is infused with ceramide



- Ceramide is a natural component of the human skin
- Ceramide links the cells of the stratum corneum together to form a waterproof, protective barrier
- Ceramide helps prevent water loss that can lead to skin damage and dryness