

HOYER® SLINGS

Quickfit Deluxe Sling



SAFE WORKING LOAD:
500 lb / 227 kg



DEPENDENCY:
Complete Dependence
FIM 1-2



ATTACHMENT:
Loops
6-Point Spreader Bar



MATERIAL OPTIONS:
Padded Polyester
Mesh

The Quickfit Deluxe sling is suitable for 95% of patients. It is available with head support and can be used for amputees following a careful assessment.

AREA	X SMALL (inch/mm)	SMALL (inch/mm)	MEDIUM (inch/mm)	LARGE (inch/mm)	X LARGE (inch/mm)
A	21.5 / 546	30.3 / 770	35 / 890	41 / 1041	52 / 1320
B	18.75 / 476	24.8 / 630	30.3 / 770	33.5 / 850	41.3 / 1049
C	27.75 / 705	36.2 / 920	44.1 / 1120	47.5 / 1206	59.3 / 1506

* XS only applies to polyester, XL only applies to mesh and polyester



Long Seat Sling



SAFE WORKING LOAD:
500 lb / 227 kg



DEPENDENCY:
Complete Dependence
FIM 1-2



ATTACHMENT:
Loops
6-Point Spreader Bar



MATERIAL OPTIONS:
Polyester
Mesh

The Long Seat sling is designed for amputees, but it is also suitable for other patients who have to remain in the sling for an extended period, i.e. during bed making when bathing.

AREA	SMALL (inch/mm)	MEDIUM (inch/mm)	LARGE (inch/mm)	X LARGE (inch/mm)
A	23.2 / 590	27.1 / 690	28.7 / 730	30.3 / 768
B	40.2 / 1020	42.0 / 1067	46.5 / 1181	49.8 / 1263

* XL applies to polyester

