

Passy-Muir Valves

All Passy-Muir® Valves

- Have the patented closed position "no leak" design
- Are designed to fit the universal 15mm hub of disposable tracheostomy tubes including neonatal and pediatric tubes
- Can be used interchangeably by tracheostomized and ventilator dependent patients with the exception of the PMV[®] 2020



PMV® 007 (Aqua Color™) Designed for ventilator application with disposable tubing



PMV® 2000 (clear) Low profile design with lower resistance, recommended for outpatient use



PMV® 2001 (Purple Color™) Low profile design with lower resistance, recommended for inpatient use



PMV® 005 (white) David's Original Design

Passy-Muir Valve Accessories

Passy-Muir® Valve Oxygen Adapter PMA® 2000

The PMA® 2000 Oxygen Adapter snaps onto the PMV® 2000 and PMV® 2001 and allows for easy delivery of low flow supplemental oxygen and humidity.



PMV® Secure-It®
The PMV® Secure-It®
strap is designed for
use with the lower
resistance 2000
series valves
and attaches

series valves and attaches the valve to the trach tie to prevent loss.



Passy-Muir Valves for Metal Tubes

PMV® 2020 (clear)

The PMV® 2020 (clear) with the PMA® 2020-S adapter are for use with metal Jackson Improved Tracheostomy Tubes (Sizes 4, 5, and 6).

PMV® 2020 .

(clear)

PMA® 2020-S



The Jackson Improved Tracheostomy Tube with Permanent 15mm Adapter can accommodate all Passy-Muir® Valves except the PMV® 2020.



Ventilator Adapters

Passy Muir Adapter Part # PMV-AD1522™ 15mm x 22mm Adapter www.passy-muir.com 800-634-5397



Passy Muir Adapter Part # PMV-AD22™ 22 mm Silicone Adapter www.passy-muir.com 800-634-5397



AirLife® Omni-Flex™ Adapter



Dual-Axis Swivel Adapter





www.passy-muir.com