

Welch Allyn Blood Pressure Cuffs – New FlexiPort Part Numbers

4/24/2008



SPOT



SPOT LXi



VSM 300



Atlas



Hand



Wall / Mobile



Pocket

Use the following part numbers for complete, ready-to-use, individual replacement cuff(s) for the devices above.

Small Infant	n/a	n/a	REUSE-06-1SC	n/a	REUSE-06-1TP	REUSE-06-2TP	REUSE-06-2BV
Infant	n/a	REUSE-07-2MQ	REUSE-07-1SC	n/a	REUSE-07-1TP	REUSE-07-2TP	REUSE-07-2BV
Small Child	REUSE-08-1SC	REUSE-08-2MQ	REUSE-08-1SC	REUSE-08-1SC	REUSE-08-1TP	REUSE-08-2TP	REUSE-08-2BV
Child	REUSE-09-1SC	REUSE-09-2MQ	REUSE-09-1SC	REUSE-09-1SC	REUSE-09-1TP	REUSE-09-2TP	REUSE-09-2BV
Small Adult	REUSE-10-1SC	REUSE-10-2MQ	REUSE-10-1SC	REUSE-10-1SC	REUSE-10-1TP	REUSE-10-2TP	REUSE-10-2BV
Adult	REUSE-11-1SC	REUSE-11-2MQ	REUSE-11-1SC	REUSE-11-1SC	REUSE-11-1TP	REUSE-11-2TP	REUSE-11-2BV
Adult Long	REUSE-11i-1SC	REUSE-11i-2MQ	REUSE-11i-1SC	REUSE-11i-1SC	REUSE-11i-1TP	REUSE-11i-2TP	REUSE-11i-2BV
Large Adult	REUSE-12-1SC	REUSE-12-2MQ	REUSE-12-1SC	REUSE-12-1SC	REUSE-12-1TP	REUSE-12-2TP	REUSE-12-2BV
Large Adult Long	REUSE-12L-1SC	REUSE-12L-2MQ	REUSE-12L-1SC	REUSE-12L-1SC	REUSE-12L-1TP	REUSE-12L-2TP	REUSE-12L-2BV
Thigh	REUSE-13-1SC	REUSE-13-2MQ	REUSE-13-1SC	REUSE-13-1SC	REUSE-13-1TP	REUSE-13-2TP	REUSE-13-2BV

If customer wishes to CONVERT TO FLEXIPORT, select tubing from Section #1 and all desired cuff(s) from Section #2

Section #1: Select the appropriate TUBE type below. Each pack includes 10 tubes with connectors...enough for 10 devices

Tube	1-SC	2-MQ	1-SC	1-SC	1-TPS	2-TPS	2-BVL
Section #2: Choose all desired cuff sizes. These are 'raw' cuffs. Notice that the part numbers are the same for all devices. That's why it's so easy !!							
Small Infant	REUSE-06	REUSE-06	REUSE-06	REUSE-06	REUSE-06	REUSE-06	REUSE-06
Infant	REUSE-07	REUSE-07	REUSE-07	REUSE-07	REUSE-07	REUSE-07	REUSE-07
Small Child	REUSE-08	REUSE-08	REUSE-08	REUSE-08	REUSE-08	REUSE-08	REUSE-08
Child	REUSE-09	REUSE-09	REUSE-09	REUSE-09	REUSE-09	REUSE-09	REUSE-09
Small Adult	REUSE-10	REUSE-10	REUSE-10	REUSE-10	REUSE-10	REUSE-10	REUSE-10
Adult	REUSE-11	REUSE-11	REUSE-11	REUSE-11	REUSE-11	REUSE-11	REUSE-11
Adult Long	REUSE-11L	REUSE-11L	REUSE-11L	REUSE-11L	REUSE-11L	REUSE-11L	REUSE-11L
Large Adult	REUSE-12	REUSE-12	REUSE-12	REUSE-12	REUSE-12	REUSE-12	REUSE-12
Large Adult Long	REUSE-12L	REUSE-12L	REUSE-12L	REUSE-12L	REUSE-12L	REUSE-12L	REUSE-12L
Thigh	REUSE-13	REUSE-13	REUSE-13	REUSE-13	REUSE-13	REUSE-13	REUSE-13

Fully Configured Cuff



Color-Coded Primary Sizes



FlexiPort Connector



Raw Cuff – Needs Tube

