

WELCH ALLYN

FlexiPort® Blood Pressure Cuffs

THE TRUE STANDARD IN BLOOD PRESSURE CUFF STANDARDIZATION.





Why choose FlexiPort®?

The true standard in blood pressure cuff standardization.

When it comes to standardizing blood pressure cuffs, Welch Allyn FlexiPort blood pressure cuffs can help hospitals **improve financial outcomes** and **enhance patient safety** with options that help streamline workflow and inventory management, and reduce the risk of cross-contamination.

The secret to true cuff standardization is the single-point FlexiPort connection, a practical innovation that enables virtually any device to work with Welch Allyn FlexiPort blood pressure cuffs—without the hassle of traditional cuff tubes and connectors.

How can FlexiPort help you?



Enhance Patient Safety

By allowing cuffs to be used with virtually any device, FlexiPort makes it easier to find the proper-sized cuff for your patient and help prevent miscuffing.





Increase Staff Satisfaction

The single-point FlexiPort connection can improve nursing workflows, delivering a more reliable connection system and fewer complaints from hospital staff.

Biomedical Engineers



Improve the Bottom Line

FlexiPort standardization can help reduce excess inventory costs by up to \$140k, lowering cuff inventory levels as much as 90%.*

Purchasing



Reduce KISK

FlexiPort technology offers hospitals options that facilitate single-patient-use to help reduce the risk of cross-contamination.



The most comprehensive cuff options to match your needs.

















True standardization—time after time.

Affordable and easy to use, Welch Allyn FlexiPort Reusable Blood Pressure Cuffs offer all the advantages of FlexiPort technology—including simplified connectivity—in a dependable, durable design.



- Tested to withstand demanding clinical environments
- Three-year warranty



 Rotatable port reduces stress to cuff tubing and port, improving patient comfort and cuff durability





Simplify and Save with Standardization

FlexiPort®: The Only True Standardization Solution

Other current cuff standardization methods typically involve a confusing array of cuff configurations, Y-tubes, clamps, and connectors. This translates into large inventory levels, confusion for the caregiver, and potential device performance issues due to incorrect or leaky connections.

With FlexiPort, Welch Allyn offers a true cuff standardization solution for your facility, allowing you to use one cuff on most electronic and manual blood pressure devices. FlexiPort cuffs not only streamline workflows by making it easier for caregivers to find and change the right cuff, but also reduce part numbers by up to 90 percent, simplifying inventory and ordering.

Why Standardize with FlexiPort?

- Streamline SKUs and inventory
- Improve consistency across your facility
- Helps clinicians simplify workflows
- More cuff options than other providers
- Options that help reduce environmental impact







ORDERING INFORMATION

FlexiPort Information: Fully Configured Cuffs

		One-Tube			Two-Tube				
					n columns		[i]oo		7
Size (20 cuffs/case)	Range	Screw	Tri-Purpose	Locking	Bayonet	Screw	Locking	Tri-Purpose & Empty Tube	Inflation Sys.
Small Infant (6)	7-10 cm	SOFT-06-1SC	SOFT-06-1TP	SOFT-06-1MQ	SOFT-06-1HP	n/a	n/a	n/a	n/a
Infant (7)	9-13 cm	SOFT-07-1SC	SOFT-07-1TP	SOFT-07-1MQ	SOFT-07-1HP	SOFT-07-2SC	SOFT-07-2MQ	SOFT-07-2TP	SOFT-07-2BV
Small Child (8)	12-16 cm	SOFT-08-1SC	SOFT-08-1TP	SOFT-08-1MQ	SOFT-08-1HP	SOFT-08-2SC	SOFT-08-2MQ	SOFT-08-2TP	SOFT-08-2BV
Child (9)	15-21 cm	SOFT-09-1SC	SOFT-09-1TP	SOFT-09-1MQ	SOFT-09-1HP	SOFT-09-2SC	SOFT-09-2MQ	SOFT-09-2TP	SOFT-09-2BV
Small Adult (10)	20-26 cm	SOFT-10-1SC	SOFT-10-1TP	SOFT-10-1MQ	SOFT-10-1HP	SOFT-10-2SC	SOFT-10-2MQ	SOFT-10-2TP	SOFT-10-2BV
Adult (11)	25-34 cm	SOFT-11-1SC	SOFT-11-1TP	SOFT-11-1MQ	SOFT-11-1HP	SOFT-11-2SC	SOFT-11-2MQ	SOFT-11-2TP	SOFT-11-2BV
Adult Long (11L)	25-34 cm	SOFT-11L-1SC	SOFT-11L-1TP	SOFT-11L-1MQ	SOFT-11L-1HP	SOFT-11L-2SC	SOFT-11L-2MQ	SOFT-11L-2TP	SOFT-11L-2BV
Large Adult (12)	32-43 cm	S0FT-12-1SC	SOFT-12-1TP	SOFT-12-1MQ	SOFT-12-1HP	SOFT-12-2SC	SOFT-12-2MQ	SOFT-12-2TP	SOFT-12-2BV
Large Adult Long (12L)	32-43 cm	SOFT-12L-1SC	SOFT-12L-1TP	SOFT-12L-1MQ	SOFT-12L-1HP	SOFT-12L-2SC	SOFT-12L-2MQ	SOFT-12L-2TP	SOFT-12L-2BV
Thigh (13)	40-55 cm	S0FT-13-1SC	SOFT-13-1TP	S0FT-13-1MQ	SOFT-13-1HP	SOFT-13-2SC	S0FT-13-2MQ	SOFT-13-2TP	SOFT-13-2BV

For FlexiPort Reusable cuffs, replace 'SOFT' with 'REUSE' in part number. For Vinyl cuffs, replace 'SOFT' with 'VINYL' in part number.

Neonatal Soft Disposable Cuffs with Male Luer Slips

Size	Range	One-Tube (10 pack)	Two-Tube (40 pack)
NEO 1	8.2-15.3 cm	39194	39114
NEO 2	7-13 cm	39193	39113
NEO 3	5.8-10.8 cm	39192	39112
NEO 4	4.6-8.5 cm	39191	39111
NEO 5	3.4-6.3 cm	39190	39110

For more information, contact your local Welch Allyn representative or visit us online at **www.welchallyn.com/cuffs**.

