

# microdot® BLEACH WIPE



**Ready to use, EPA registered  
hospital use disinfectant bleach wipe**

**SPORICIDAL • BACTERICIDAL • FUNGICIDAL • \*VIRUCIDAL**  
**Disinfects • Cleans • Deodorizes • Kills 99.99% of Germs**

## ONE STEP CLEANER AND DISINFECTANT THAT KILLS:

- C DIFF Spores in ONE STEP in 3 minutes
- HBV in 30 seconds
- HCV in 30 seconds (preclean prior to disinfecting surface)

## KILLS THE FOLLOWING ORGANISMS IN 5 MINUTES:

- Methicillin-Resistant *Staphylococcus aureus* (MRSA)
- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- Vancomycin-Resistant *Enterococcus faecalis* (VRE)
- Human Immunodeficiency Virus Type 1 (HIV-1)
- *Trichophyton mentagrophytes*

Premoistened wipes with a 1:10 dilution of bleach following the guidelines of the CDC, OSHA, CMS, SHEA-IDSA and APIC

CAN BE USED ON BLOOD GLUCOSE METERS.

† Disinfects hard, non-porous healthcare surfaces, patient care equipment and point-of-care equipment.

† This product is for the exterior surfaces of Blood Glucose Meters and is not for use on the needle and/or the monitors.