



Monoject™ Phlebotomy Sharps Containers

Engineering features include:

- Fit in most common caddies and trays
- Available holders and brackets
- A variety of container sizes and openings allow for disposal of other sharps such as lancets, butterfly tubing and small syringes
- Containers lock for final disposal

Order code	Description	Container dimensions H x D x W in.	Case qty.
8900SA	1 quart, red	6.25 x 4.25 x 4.25	100
1522SA	2.2 quart, red	7.25 x 4.75 x 6.5	60
8900MW	Renewable, 1 quart, opaque red	6.25 x 4.25 x 4.5	100



Monoject™ Large Volume Sharps Containers

Engineering features include:

- Sliding, hinged and sealing gasket lid
- Adjustable sliding and hinged lids accommodate large sharps
- Lids can provide temporary closure when container is not in use to reduce exposure to contents
- Carrying handles and carts with wheels are available for easy transport and lifting of heavy, full containers
- Optional gasket lids provide leak-resistant containment of contents
- Easy, simple assembly
- Containers lock for final disposal
- Nestable containers save valuable storage space
- Sliding gasket lids have a round 6 in. opening to accommodate up to 1,000 mL bottles



Monoject™ Safety Large Volume Sharps Containers

Order code	Description	Container dimensions H x D x W in.	Lid opening dimensions L x W in.	Case qty.
31143665	10 gallon, red, split lid	21.5 x 12 x 15.5	Rectangular, 8.5 x 1.75	6

31143665

Reminder: Containers may vary in color, but are most generally indicated by colors in construction.

Biohazardous