Therapy Putty

Prepacked let us do the work!



5 lbs. of putty packed into 40 2-oz. containers

■ Includes 40 free containers

Putty Packs save you precious time and money. We've taken 5 lbs. of putty and packaged it into 2-oz. containers at no additional cost. We give you the 40 containers free and save you the time and hassle of packaging your own putty. Not made with natural rubber latex.

| #2970001342 | TAN, XX-SOFT | 40/PACK | |
|-------------|----------------|---------|--|
| #2970001344 | YELLOW, X-SOFT | 40/PACK | |
| #2970001346 | RED, SOFT | 40/PACK | |
| #2970001348 | GREEN, MEDIUM | 40/PACK | |
| #2970001350 | BLUE, FIRM | 40/PACK | |

• 55 California residents, see below.

Visit

AliMed.com/putty-packs

for more details

Standard putty-for less!

■ Guaranteed to maintain consistency for an entire year



Color-Coded Therapy Putty comes in resistances consistent with industry standards. Silicone-based putty. Guaranteed to maintain its consistency for one year from the date you receive it. Not made with natural rubber latex.







| COLOR/RESISTANCE | 2 OZ. | 3 OZ. | 4 OZ. | 1 LB. | 5 LBS. |
|------------------|-------------|-------------|-------------|-------------|-------------|
| TAN, XX-SOFT | #2970001341 | #2970001362 | #2970001367 | #2970001372 | #2970001377 |
| YELLOW, X-SOFT | #2970001343 | #2970001363 | #2970001368 | #2970001373 | #2970001378 |
| RED, SOFT | #2970001345 | #2970001364 | #2970001369 | #2970001374 | #2970001379 |
| GREEN, MED | #2970001347 | #2970001365 | #2970001370 | #2970001375 | #2970001380 |
| BLUE, FIRM | #2970001349 | #2970001366 | #2970001371 | #2970001376 | #2970001381 |
| | EACH | EACH | EACH | EACH | EACH |

■ • SS California residents, see below.

PERSONAL PUTTY KIT

#2970003552 ONE OF EACH: 2-oz PACK OF YELLOW, RED, GREEN, AND BLUE 4 PACKS/KI

• 65 California residents, see below.

Containers



■ Available in three sizes

Putty Containers help maintain the consistency of putty. Available in three sizes.

| #2970002754 | 2-oz | 10/PACK |
|-------------|------|---------|
| #2970002752 | 4-oz | 10/PACK |
| #2970002734 | 6-oz | 10/PACK |

65 California residents, see below.

■ Simulate everyday activities

Puttycise® Tools can be pushed, pulled, or turned through putty to produce specified exercises. Changing the putty resistance increases or decreases difficulty level. Use to simulate functional activity and/or teach joint protection techniques. Consists of five tools that all function by inserting tool base into resistive putty. Tools include Knob Turn, Peg Turn, Key Turn, Cap Turn, and L-Bar (simulates pulling a vacuum, sweeping, using a hammer, etc.). Tools work best in 1 lb. of putty.

5-Piece Set includes one of each tool with a carry bag and set of instruction cards. Putty not included.

| #2970000935 | 5-PIECE SET w/CARDS AND BAG | SET |
|-------------|-----------------------------|------|
| #2970000936 | KNOB TURN | EACH |
| #2970000938 | PEG TURN | EACH |
| #2970000939 | KEY TURN | EACH |
| #2970000940 | CAPTURN | EACH |
| #2970000937 | L-BAR | EACH |
| | | |

65 California residents, see below.







Knob Turn Peg Turn Key





Scan to shop more Therapy Supplies online!

AliMed.com/exercise

