

Filac™ 3000 EZ Electronic Thermometer

- Icon based user interface on backlit LCD display
- Push button alternates Fahrenheit/Celsius
- Interchangeable, color-coded isolation chambers with matching probes
- Auto on™, “Auto off” feature
- Last temperature recall
- Utilizes item 202020, 502000, or 500500 disposable probe cover

Order Information

Item code	Description	Ship case
504000	Filac™ 3000 EZ Oral/Axillary Complete	1
504005	Filac™ 3000 EZ Rectal Complete	1
504027	Filac™ 3000 EZ Oral/Axillary Complete – 9ft	1
504037	Filac™ 3000 Rectal Complete – 9 ft	1
500026	Filac™ 3000 Oral Probe – 4 ft	1
500027	Filac™ 3000 Oral Probe – 9ft	1
500028	Filac™ 3000 Oral Isolation Chamber	1
500036	Filac™ 3000 Rectal Probe – 4ft	1
500037	Filac™ 3000 Rectal Probe – 9 ft	1
500038	Filac™ 3000 Rectal Isolation Chamber	1
500099	Filac™ 3000 Calibration Plug	1

Filac 3000 EZ hardware is warranted for three years from shipment for manufacturer’s defect only.

PRODUCT SPECIFICATIONS

Accuracy

- Oral Quick Mode: $\pm 0.3^{\circ}\text{C}$
- Oral Standard Mode: $\pm 0.1^{\circ}\text{C}$
- Axillary: $\pm 0.1^{\circ}\text{C}$
- Rectal: $\pm 0.1^{\circ}\text{C}$

Measurement response time

- Oral Quick Mode
 - 3-5 seconds for non-fever temperatures
 - 8-10 seconds for fever temperatures
- Axillary: 8–12 seconds
- Rectal: 10–14 seconds
- Direct Mode: 60–120 seconds

Temperature measurement range

- $30\text{--}43^{\circ}\text{C}$ ($86^{\circ}\text{F}\text{--}109^{\circ}\text{F}$)

Battery life

- 6000 temperatures (4 AA batteries)



504000



504005



500036