

## PolyMem® and PolyMem Silver® Non-Adhesive Dressings

### PolyMem Non-Adhesive

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
5022	1.8" x 1.8" (4.7cm x 4.7cm)	100	5 bx / 20 ea	A6209
5033	3" x 3" (8cm x 8cm)	60	4 bx / 15 ea	A6209
5044	4" x 4" (10cm x 10cm)	60	4 bx / 15 ea	A6209
5055	5" x 5" (13cm x 13cm)	30	2 bx / 15 ea	A6210
5077	6.5" x 7.5" (17cm x 19cm)	30	2 bx / 15 ea	A6210
5124	4" x 12.5" (10cm x 32cm)	12	1 bx / 12 ea	A6210
5244	4" x 24" (10cm x 61cm) Roll	8	2 bx / 4 ea	A6211
5824	8" x 24" (20cm x 60cm) Roll	2	1 bx / 2 ea	A6211

### PolyMem Silver Non-Adhesive

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
1022	1.8" x 1.8" (4.7cm x 4.7cm)	100	5 bx / 20 ea	A6209
1044	4.25" x 4.25" (10.8cm x 10.8cm)	30	2 bx / 15 ea	A6210
1077	6.5" x 7.5" (17cm x 19cm)	15	1 bx / 15 ea	A6211
1124	4.25" x 12.5" (10.8cm x 32cm)	12	1 bx / 12 ea	A6211

# PolyMem®

# PolyMem® SILVER

The optimal choice when cutting dressings to size and tapeless attachment or roll taping, is your preferred method of application. Use as a primary, secondary or combined primary and secondary dressing.



## PolyMem MAX and PolyMem MAX Silver Non-Adhesive Dressings

### PolyMem MAX


REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
5035	3" x 3" (7.6cm x 7.6cm)	40	4 bx / 10 ea	A6209
5045	4.5" x 4.5" (11cm x 11cm)	20	2 bx / 10 ea	A6210
5088	8" x 8" (20cm x 20cm)	10	2 bx / 5 ea	A6211

### PolyMem MAX Silver

1045	4" x 4" (10cm x 10cm)	16	2 bx / 8 ea	A6209
1088	8" x 8" (20cm x 20cm)	10	2 bx / 5 ea	A6211

# PolyMem® MAX®

This configuration is a thicker PolyMem pad with high MVTR film for wounds with heavier drainage or when longer wear time is desired. Use as a primary, secondary or combined primary and secondary dressing.



## PolyMem WIC, PolyMem WIC Silver, and PolyMem WIC Silver Rope Wound Fillers

### PolyMem WIC

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
5733	3" x 3" (8cm x 8cm) 4 grams	40	4 bx / 10 ea	A6215
5712	3" x 12" (8cm x 30cm) 16 grams	12	1 bx / 12 ea	A6215
5788	8" x 8" (20cm x 20cm) 28 grams	10	2 bx / 5 ea	A6215

### PolyMem WIC Silver

1331	1" x 3" (2.5cm x 7.6cm) 1.3 grams	70	5 bx / 14 ea	A6215
1333	3" x 3" (8cm x 8cm) 4 grams	20	2 bx / 10 ea	A6215
1788	8" x 8" (20cm x 20cm) 28 grams	10	2 bx / 5 ea	A6215

### PolyMem WIC Silver Rope

1814	0.4" x 14" (1cm x 35cm) 3 grams	12	2 bx / 6 ea	A6215
------	---------------------------------	----	-------------	-------

# PolyMem® WIC®

With no film backing, WIC expands within the wound cavity to fill dead space, keeping the wound bed clean and moisturized while absorbing excess fluid. WIC Rope has a reinforced mesh to ensure one-piece removal.



## PolyMem Finger/Toe Dressings

### Finger / Toe

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
4401	#1 SM; circ. 1.8" - 2.2" (4.67cm - 5.7cm)	30	5 bx / 6 ea	A6209
4402	#2 MED; circ. 2.2" - 2.6" (5.7cm - 6.72cm)	30	5 bx / 6 ea	A6209
4403	#3 LG; circ. 2.6" - 3.0" (6.72cm - 7.74cm)	30	5 bx / 6 ea	A6209
4404	#4 XL; circ. 3.0" - 3.4" (7.74cm - 8.76cm)	30	5 bx / 6 ea	A6209
4405	#5 XXL; circ. 3.4" - 3.8" (8.76cm - 9.78cm)	30	5 bx / 6 ea	A6210

### Finger / Toe Silver

REF	DIMENSIONS	PER CASE	PACKAGING	HCPCS Codes
1401	#1 SM; circ. 1.8" - 2.2" (4.67cm - 5.7cm)	30	5 bx / 6 ea	A6209
1402	#2 MED; circ. 2.2" - 2.6" (5.7cm - 6.72cm)	30	5 bx / 6 ea	A6209
1403	#3 LG; circ. 2.6" - 3.0" (6.72cm - 7.74cm)	30	5 bx / 6 ea	A6209
1404	#4 XL; circ. 3.0" - 3.4" (7.74cm - 8.76cm)	30	5 bx / 6 ea	A6209
1405	#5 XXL; circ. 3.4" - 3.8" (8.76cm - 9.78cm)	30	5 bx / 6 ea	A6210

This configuration is specially designed for use on fingers and toes. Appropriate for use as a combined primary and secondary dressing or as a secondary dressing.

