


CATHETER SECUREMENT

The Revolution™ Catheter Securement Device saves time, provides convenience and enhances compliance with safety regulations for securing catheters with a needleless procedure. The StayFIX® Fixation Device design secures catheters in place without the use of sutures, reducing the exposure to sharps and sharps-related accidents.



REVOLUTION™

Catheter Securement Device

- Secures non-vascular catheters from 6F - 24F
- Hydrocolloid is absorbent around puncture site and conforms well to body anatomy
- Provides safe and secure fixation for up to 7 days
- Included shower cap covers insertion site for added protection while showering

Catalog Number	Description
CFG624	6F–24F

Packaged 25 per box

StayFIX®

Fixation Device

- Secures non-vascular catheters from 5F–22F
- Provides safe and secure fixation for up to 7 days

Catalog Number	Description
680ME*	5F–14F
685ME**	12F–22F

* Packaged 25 per box | ** Packaged 20 per box

