

3M™ Tegaderm™ PICC/CVC Dressing Change Kits

Product Information

3M | Tegaderm™ PICC/CVC

REF 77-2100K01

Dressing Change Kit

Contents:

- 1 Tegaderm™ CHG w/Device
- 1 Cavilon™ No Sting Barrier
- 1 Anti: Chloraprep® 3ml w/insert
- 1 Towel: 13X18 2ply Folder
- 1 Mask w/Earloop Folder

Country
of Origin
US
US
US
CN
CN

- 2 Gauze: 4X4 4-Ply NW
- 1 Tape: Measure 24 in
- 1 Anti: Alcohol Swab 3
- 1 Glove: Nitrile Med Pf Blu Wal

Country
of Origin
CN
US
US
LK

US – United States, CN – China, LK-Sri Lanka

CAUTION: Store at room temperature.

Warning! Extremely Flammable.

Caution: Federal law restricts the device to sale by or on the order of a licensed health care professional.

CAUTION: Sterilant does not penetrate glass ampules or foil packets.

ChorlaPrep® is a registered trademark of CareFusion

ChorlaPrep® - www.carefusion.com



Promoting

Consistent Care

In one convenient package, Tegaderm™ PICC/CVC Dressing Change Kits include 3M's latest products along with all necessary components for securing and protecting catheters, and maintaining skin integrity.

To meet diverse clinical needs, six varieties of dressing change kits are available. Each kit features a 3M™ PICC/CVC Securement Device and a securement dressing, a choice of two antimicrobial CHG dressings, and additional components including gloves, masks, CHG preps, gauze, and more.

Tegaderm™
Brand

3M

3M™ PICC/CVC Securement System

Each kit contains a version of the 3M™ PICC/CVC Securement System, a sutureless securement system that balances secure adhesion with gentleness to skin. The system is made up of the 3M™ PICC/CVC Securement Device and one of two Tegaderm™ securement dressings: Tegaderm™ I.V. Advanced Securement Dressing or Tegaderm™ CHG Chlorhexidine Gluconate I.V. Securement Dressing.

- Designed to be worn for up to 7 days
- Skin-friendly materials improve comfort

- Removes cleanly without alcohol, and without causing patients undue pain or stress
- Accommodates the majority of PICC and CVC catheters up to and including 12 French
- Device uses a silicone adhesive, allowing for strong hold yet gentle removal
- Integrated antimicrobial CHG gel pad of Tegaderm™ CHG provides immediate and continuous protection

3 kits will feature:

- 3M™ PICC/CVC Securement Device
- Tegaderm™ CHG Securement Dressing
- 3M™ Caviion™ No Sting Barrier Film
- ChloraPrep® solution

Catalog Number	Included Items	Kits/Case
77-2100K01	• 3M™ PICC/CVC Securement Device (SKU # 2100) • Tegaderm™ CHG Dressing (SKU # 1877) • Caviion™ No Sting Barrier Film (SKU # 3343K) • Gauze: 4 inches x 4 inches (qty 2) • Gloves – Powder-Free, Nitrile, Blue (qty 1)	25
77-2100K03	• 3M™ PICC/CVC Securement Device (SKU # 2100) • Tegaderm™ CHG Dressing (SKU # 1877) • Caviion™ No Sting Barrier Film (SKU # 3343K) • Gauze: 4 inches x 4 inches (qty 2) • Gloves – Powder-Free, Nitrile, Blue (qty 1)	25
79-2100K01 (features largest Tegaderm™ CHG Dressing)	• 3M™ PICC/CVC Securement Device (SKU # 2100) • Tegaderm™ CHG Dressing (SKU # 1879) • Caviion™ No Sting Barrier Film (SKU # 3343K) • Gauze: 4 inches x 4 inches (qty 2) • Gloves – Powder-Free, Nitrile, Blue (qty 1)	25

3 kits will feature:

- 3M™ PICC/CVC Securement Device
- Tegaderm™ I.V. Advanced Dressing
- 3M™ Caviion™ No Sting Barrier Film
- ChloraPrep® solution
- BioPatch® Protective Disk

Catalog Number	Included Items	Kits/Case
37-2100K01	• 3M™ PICC/CVC Securement Device (SKU # 2100) • Tegaderm™ I.V. Advanced Dressing (SKU # 1837) • Caviion™ No Sting Barrier Film (SKU # 3343K) • Gauze: 4 inches x 4 inches (qty 2) • BioPatch® Protective Disk w/CHG (qty 1)	25
37-2100K03	• 3M™ PICC/CVC Securement Device (SKU # 2100) • Tegaderm™ I.V. Advanced Dressing (SKU # 1837) • Caviion™ No Sting Barrier Film (SKU # 3343K) • Gauze: 4 inches x 4 inches (qty 2) • BioPatch® Protective Disk w/CHG (qty 1)	25
39-2100K01 (features largest Tegaderm™ I.V. Advanced Dressing)	• 3M™ PICC/CVC Securement Device (SKU # 2100) • Tegaderm™ I.V. Advanced Dressing (SKU # 1839) • Caviion™ No Sting Barrier Film (SKU # 3343K) • Gauze: 4 inches x 4 inches (qty 2) • BioPatch® Protective Disk w/CHG (qty 1)	25



Critical & Chronic Care Solutions Division

3M Health Care
2510 Conway Avenue
St. Paul, MN 55144
USA
1-800-228-3957
www.3M.com/C3S

Please recycle. Printed in U.S.A.
© 3M 2015. All rights reserved.
3M, Caviion and Tegaderm are trademarks of 3M.
Used under license in Canada. BioPatch and
ChlorPrep are not trademarks of 3M.
70-2011-5646-3

Tegaderm™
Brand