



Monoject™ Transportable Sharps Containers

Engineering features include:

- Simple, one-handed operation opens and closes the one-quart disposal opening
- Convenient size fits into carry bags for nursing, EMS and fire/rescue
- Large disposal opening easily accommodates small syringes and other sharps waste
- Sharps shuttle snap closure for one-handed opening and closing

Order code	Description	Container dimensions H x D x W in.	Case qty.
8303SA	1 quart, red	8.75 x 2.5 x 4.5	20



8303SA



Reminder: Containers may vary in color, but are most generally indicated by colors in construction.

■ Biohazardous