

## Alginates



### Calcium Alginate Dressings – Kendall™

- Absorbs exudate from wound and creates a protective cushioning gel and optimal moist healing environment
- Protective gel binds exudate and helps prevent skin maceration
- Can be used on infected wounds
- Absorbs approximately 20 times the dressing's weight
- Unique composition and weaving pattern help provide a simplified one piece removal
- Not made with natural rubber latex

#### Calcium Alginate Dressings

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
9232	2 in. x 2 in.		Sterile	1 ea/po, 10 po/bx, 10 bx/cs	100
9233-	4 in. x 4 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
9236	4 in. x 4 in.	Plus	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
9240	4 in. x 5 ½ in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
9238	4 in. x 8 in.		Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
9239	6 in. x 10 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
9242	12 in. x 24 in.		Sterile	1 ea/po, 5 po/cs	5
9231-	12 in.	Rope	Sterile	1 ea/po, 5 po/bx, 4 bx/cs	20
9243	24 in.	Rope	Sterile	1 ea/po, 5 po/bx, 4 bx/cs	20
9244	36 in.	Rope	Sterile	1 ea/po, 5 po/bx, 4 bx/cs	20

#### Zinc Calcium Alginate Dressings

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
9354-	2 in. x 2 in.	with Zinc	Sterile	1 ea/po, 10 po/bx, 10 bx/cs	100
9355-	4 in. x 4 in.	with Zinc	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
9356-	4 in. x 8 in.	with Zinc	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
9351-	12 in.	Rope with Zinc	Sterile	1 ea/po, 5 po/bx, 4 bx/cs	20



### Reinforced Gelling Fibers (Alginate) – Cardinal Health™

- Non-woven calcium/sodium alginate consisting of a high M alginate in a cohesive fiber construction
- Our reinforced gelling fibers assist in maintaining a moist environment for optimal wound healing and the formation of granulation tissues
- Highly absorbent gelling fibers wick and retain exudate away from the wound bed
- Cohesive fiber construction facilitates easy intact/one piece atraumatic removal
- The silver ions in our reinforced gelling Fibers+Ag protect the dressing from a broad spectrum of microorganisms over a period of up to 14 days
- Not made with natural rubber latex

#### Reinforced Gelling Fibers

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
RGFB22	2 in. x 2 in.		Sterile	1 ea/po, 10 po/bx, 10 bx/cs	100
RGFB45	4 in. x 4 ¾ in.		Sterile	1 ea/po, 10 po/bx, 10 bx/cs	100
RGFB66	6 in. x 6 in.		Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
RGFBROPE	18 in.	Rope	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50

#### Reinforced Gelling Fibers +AG

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
RGFB22AG	2 in. x 2 in.	+AG	Sterile	1 ea/po, 10 po/bx, 10 bx/cs	100
RGFB45AG	4 in. x 4 ¾ in.	+AG	Sterile	1 ea/po, 10 po/bx, 10 bx/cs	100
RGFB66AG	6 in. x 6 in.	+AG	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
RGFB812AG	8 in. x 12 in.	+AG	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
RGFBROPEAG	18 in.	Rope +AG	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50