

Ortho Upper > Slings & Shoulder Supports

Rolyan® Hemi Arm Slings

Effective unilateral anti-subluxation support



Rolyan® Hemi Arm Slings.

- Positions the head of the humerus in the glenoid fossa and the scapula in external rotation.
- Enhances normal gait pattern.
- Elastic straps are dynamic, allowing a patient to actively move, return muscles and perform self-ranging exercises. The supports enable hemiplegics to use their affected arm as an assist during functional tasks.
- Easy on/off and can be worn under clothing.
- Effective following stroke, soft tissue injury or ligament strain.

A. Original Rolyan® Hemi Arm Sling. Cuff system of nylon with polyester foam laminate inside cuff. White with blue lining.

B. Rolyan® Hemi Arm Sling II. Approximately 1 1/2" (3.8 cm) longer than the original model. Nylon and polyester/elastic blend. Provides increased support for heavier patients as well as patients with increasing active motion.

C. Rolyan® Hemi Arm Sling Cuff Adapter. Fits all sizes, left and right. Positions the hand above heart to treat or prevent edema. Attaches to a D-ring on the unaffected side. Nylon with polyester foam laminate.

A653CA Fits Original Hemi Arm Sling & Custom Hemi Arm Sling **\$17.00**

A65311 Fits Hemi Arm Sling II **\$17.00**

A. Original		B. Hemi II		Size	Acromion / Chest
Left	Right	Left	Right		
A653LS	A653RS	A6531	A6532	Small	14"-18" (36-46 cm)
A653LM	A653RM	A6533	A6534	*Medium Reg.	18"-24" (46-61 cm)/*Chest Circ. 30"-35" (76-89 cm)
A653LML	A653RML	A6535	A6536	*Medium Long	18"-24" (46-61 cm)/*Chest Circ. 35"-40" (89-102 cm)
A653LL	A653RL	A6537	A6538	Large	24"-28" (61-71 cm)
A653LXL	A653RXL	A6539	A65310	X-Large	28"-34" (71-86 cm)
\$55.00		\$58.54			

Sizing Guidelines



To fit, measure the distance from acromion to acromion posteriorly across the midpoint of the scapula.



To determine if patient is medium or medium long, measure the chest circumference at the axilla.

***Chest measurement only needed for size medium patients.**

Do you have patients who are hard to fit? Order the cuffs and straps separately as shown below.

Rolyan® Custom Hemi Arm Sling. White durafoam cuff is lined with Scotchguard to prevent the cuff from migrating. Latex free.

Bicep Cuffs		Size	Bicep Circ.
Left	Right		
55986401	55986404	Small	10"-12" (25.4-30.5 cm)
55986402	55986405	Medium	12"-14" (30.5-35.6 cm)
55986403	55986406	Large	14"-16" (35.6-40.6 cm)
\$18.50			

Straps		Size	Acromion as shown in the top diagram of sizing guide above
Left	Right		
55986407	55986410	Small	16"-22" (40.6-55.9 cm)
55986408	55986411	Medium	22"-28" (55.9-71.1 cm)
55986409	55986412	Large	28"-34" (71.1-86.4 cm)
\$39.15			



If forearm support is needed, order A653CA Cuff Adapter as shown above in Image C.

636

Call toll-free 1.800.323.5547

fax 1.800.547.4333