

New Image Two-Piece Pouching System

Convex Skin Barriers with Integrated Floating Flange (continued)



1150x (PRE-SIZED)

- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border

Box of 5



HCPCS: A4407

Stock No	Flange Size	Barrier Opening
11503	1-3/4" (44 mm)	7/8" (22 mm)
11504	1-3/4" (44 mm)	1" (25 mm)
11505	2-1/4" (57 mm)	1-1/8" (29 mm)
11506	2-1/4" (57 mm)	1-1/4" (32 mm)

Upgraded Design Recommended replacement for 149xx series below



139xx (PRE-SIZED)

- ▶ *Flextend* (extended wear) Skin Barrier
- ▶ With Tape Border

Box of 5



Upgraded Design Features:

- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edge

HCPCS: A4407

Stock No	Flange Size	Barrier Opening
13901	1-3/4" (44 mm)	5/8" (16 mm)
13902	1-3/4" (44 mm)	3/4" (19 mm)
13903	1-3/4" (44 mm)	7/8" (22 mm)
13904	1-3/4" (44 mm)	1" (25 mm)
13905	2-1/4" (57 mm)	1-1/8" (29 mm)
13906	2-1/4" (57 mm)	1-1/4" (32 mm)
13907	2-1/4" (57 mm)	1-3/8" (35 mm)
13908	2-1/4" (57 mm)	1-1/2" (38 mm)
13909	2-3/4" (70 mm)	1-5/8" (41 mm)
13910	2-3/4" (70 mm)	1-3/4" (44 mm)
13911	2-3/4" (70 mm)	2" (51 mm)

Original Design

149xx (PRE-SIZED)

- ▶ *Flextend* (extended wear) Skin Barrier
- ▶ With Tape Border

Box of 5



HCPCS: A4407

Stock No	Flange Size	Barrier Opening
14901	1-3/4" (44 mm)	5/8" (16 mm)
14902	1-3/4" (44 mm)	3/4" (19 mm)
14903	1-3/4" (44 mm)	7/8" (22 mm)
14904	1-3/4" (44 mm)	1" (25 mm)
14905	2-1/4" (57 mm)	1-1/8" (29 mm)
14906	2-1/4" (57 mm)	1-1/4" (32 mm)
14907	2-1/4" (57 mm)	1-3/8" (35 mm)
14908	2-1/4" (57 mm)	1-1/2" (38 mm)
14909	2-3/4" (70 mm)	1-5/8" (41 mm)
14910	2-3/4" (70 mm)	1-3/4" (44 mm)
14911	2-3/4" (70 mm)	2" (51 mm)



DID YOU KNOW...

Convex Skin Barrier Options

The New Image two-piece pouching system delivers several skin barrier options to fit a variety of patient needs. Built-in convexity options provide security for recessed and flush stomas, and peristomal skin with soft muscle tone. Convex skin barriers are also available with and without tape.