

## AMPUTEE CLIP IN SITU SLING

The Amputee clip in situ sling is suitable for residents/patients with bilateral leg amputations above the knee or total hip. It is designed with the ability to remain under the patient between transfers following clinical assessment.

PRODUCT INFORMATION	
<b>ARTICLE NUMBER</b>	<b>MAA2050M</b>
Safe working load	190 kg (418 lbs)
Size	XS S M L LL XL
Morphology	
Cleaning method	
Material	Polyester woven fabric

### COMPATIBLE LIFTS

- |   |  |
|---|--|
| <b>Floor lifts</b>  | <b>Ceiling lifts</b>   |
| <ul style="list-style-type: none"> <li>Maxi Move</li> <li>Maxi Twin</li> <li>Maxi Twin compact</li> <li>Maxi 500</li> </ul> | <ul style="list-style-type: none"> <li>Maxi Sky 600</li> <li>Maxi Sky 2 Series</li> <li>Maxi Sky 2 Plus (Single mode)</li> </ul> |



## AMPUTEE CLIP BATHING IN SITU SLING

The Amputee clip bathing in situ sling is suitable for bathing and showering residents/patients with bilateral leg amputations above the knee or total hip. It is designed with the ability to remain under the patient between transfers following clinical assessment.

PRODUCT INFORMATION	
<b>ARTICLE NUMBER</b>	<b>MAA2080M</b>
Safe working load	190 kg (418 lbs)
Size	XS S M L LL XL
Morphology	
Cleaning method	
Material	Lightweight knitted polyester fabric - mesh

### COMPATIBLE LIFTS

- |   |  |
|---|--|
| <b>Floor lifts</b>  | <b>Ceiling lifts</b>   |
| <ul style="list-style-type: none"> <li>Maxi Move</li> <li>Maxi Twin</li> <li>Maxi Twin compact</li> <li>Maxi 500</li> </ul> | <ul style="list-style-type: none"> <li>Maxi Sky 600</li> <li>Maxi Sky 2 Series</li> <li>Maxi Sky 2 Plus (Single mode)</li> </ul> |



## AMPUTEE CLIP TOILET SLING

The Amputee clip toilet sling is designed with a commode hole to assist with toileting and hygiene. This sling is suitable for general transfers of residents/patients with bilateral leg amputations above the knee or total hip.

PRODUCT INFORMATION	
<b>ARTICLE NUMBER</b>	<b>MAA2091M</b>
Safe working load	272 kg (600 lbs)
Size	XS S M L XL
Morphology	
Cleaning method	
Material	Polyester knitted fabric

### COMPATIBLE LIFTS

- |   |  |
|---|--|
| <b>Floor lifts</b>  | <b>Ceiling lifts</b>   |
| <ul style="list-style-type: none"> <li>Maxi Move</li> <li>Maxi Twin</li> <li>Maxi Twin compact</li> <li>Maxi 500</li> </ul> | <ul style="list-style-type: none"> <li>Maxi Sky 600</li> <li>Maxi Sky 2 Series</li> <li>Maxi Sky 2 Plus (Single mode)</li> </ul> |



## PAEDIATRIC CLIP SLING

The Paediatric clip sling is suitable for everyday patient handling routines such as bed and chair related transfers and it is designed specifically for children with an additional buckle in the head area.

PRODUCT INFORMATION	
<b>ARTICLE NUMBER</b>	<b>MAA4020M</b>
Safe working load	125 kg (275 lbs)
Size	XXS XS S
Morphology	
Cleaning method	
Material	Polyester woven fabric

### COMPATIBLE LIFTS

- |   |  |
|---|--|
| <b>Floor lifts</b>  | <b>Ceiling lifts</b>   |
| <ul style="list-style-type: none"> <li>Maxi Move</li> <li>Maxi Twin</li> <li>Maxi Twin compact</li> <li>Maxi 500</li> </ul> | <ul style="list-style-type: none"> <li>Maxi Sky 2 Series</li> <li>Maxi Sky 2 Plus (Single mode)</li> <li>Maxi Sky 2 Series</li> <li>Maxi Sky 2 Plus (Single mode)</li> </ul> |

