# OxyMask<sup>™</sup> ETCO<sub>2</sub>

Safer care for an exceptional experience



Introducing OxyMask ETCO₂—a unique open-style oxygen mask by Southmedic. Its proprietary design allows a single mask to provide a broad range of oxygen concentrations and accurate end-tidal CO₂ monitoring.



### Improve operational efficiencies

One mask for all patients from the operating room to recovery.

- · Available with full male gas sampling line or pigtail to accommodate facility needs
- Open design allows for oral and nasal access
- Broad FiO<sub>3</sub> range (24%–65%) eliminates the need for many traditional oxygen masks
- Save caregiver time and reduce errors by reducing number of device changeouts
- Safer, more efficient oxygen delivery allows for easier patient weaning for quicker recovery



### Improve patient experience

Lightweight, open-mask design.

- Reduces claustrophobic feelings
- · Unrestricted communication with caregiver and family
- Provides easy access for fluids and oral medication delivery

## Save space in the supply room

**Before OxyMask** 



After OxyMask






One acute care facility saw a

33% SKU Reduction





# OxyMask<sup>™</sup> ETCO<sub>2</sub>

Open design mask delivers oxygen flows from 1–15 liters per minute at a range of 24%–65% FiO<sub>2</sub> with end-tidal CO<sub>2</sub> monitoring.

#### With female pigtail adapter

Item No.	Size	Universal Tubing Length	Female Pigtail Adapter Length	Pkg.
OM21258ML	Adult	7'	8"	25/bx
OP21258ML	Adult Plus	7'	8"	25/bx
OK21258ML	Pediatric	7'	8"	25/bx
OM212514ML	Adult	14'	15"	25/bx

#### With male gas sampling line

Item No.	Size	Universal Tubing Length	Male Gas Sampling Line Length	Pkg.
OM21258SLMML	Adult	7'	8'	25/bx
OP21258SLMML	Adult Plus	7'	8'	25/bx
OK21258SLMML	Pediatric	7'	8'	25/bx
OM212514SLM	Adult	14'	15'	25/bx



# Efficiency runs in the family.

Learn more about the full range of OxyMask products.

medline.com/go/oxymask



#### **Empower best practice.**

Contact your Medline Representative to learn more.



Follow us fin BLOG



#### Medline Industries, Inc.

Three Lakes Drive, Northfield, IL 60093 Medline United States | 1-800-MEDLINE (633-5463) medline.com | info@medline.com

#### Medline Canada

1-800-268-2848 | medline.ca | canada@medline.com **Medline México** 

01-800-831-0898 | medlinemexico.com | mexico@medline.com

1. Hoerr C. Whitepaper. Leveraging technology to maximize efficiencies in today's challenging respiratory care environment. Accessed April 8, 2019. Online: http://thebetteroxygenmask.com/wp-content/uploads/2016/01/Leveraging-Available-Technology-to-Maximize- Efficiencies-in-the-Provision-of-Respiratory-Carev3.pdf.