



Vista® basic Infusion Pump

Description	Product Code	Unit Per cs
Vista® basic Infusion Pump  Features include: Pump based free flow protection; Straight-line tubing; Quick start-up for continuous infusions by entering rate and volume, then press START; Program Mode and Ramp/Taper Mode; Helpful prompts to walk the user through the programming process. Vista basic Infusion Pump		
	637-202	1
Vista basic Infusion Pump with 1.2 micron filter and backcheck valve  (15 drops/mL) Universal spike, roller clamp, Supor 1.2 micron air eliminating filter, backcheck valve approx. 28 in. above distal end, standard injection site 6 in. above distal end, SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 22 mL, Length: 117 in. (approx. 297 cm)		
	352456	50
Vista basic Infusion Pump with CARESITE® Luer Access Device and backcheck valve  (15 drops/mL) Universal spike, roller clamp, backcheck valve approx. 28 in. above distal end, CARESITE® luer access device 6 in. above distal end, SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 17 mL, Length: 104 in. (approx. 264 cm)		
	352719	50
Vista basic Infusion Pump with 0.2 micron filter and backcheck valve  (15 drops/mL) Universal spike, backcheck valve, two CARESITE luer access devices approx. 80 in. and 6 in. above distal end, roller clamp, SUPOR® 0.2 micron air eliminating filter, backcheck valve and SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 23 mL, Length: 124 in. (approx. 315 cm)		
	US9923F	50
Vista basic Infusion Pump with backcheck valve  (15 drops/mL) Universal spike, roller clamp, backcheck valve approx. 28 in. above distal end, SPIN-LOCK® connector. Not made with DEHP or natural rubber latex. Priming volume: 17 mL, Length 108 in. (approx. 274 cm)		
	V9905	50
Vista basic Infusion Pump with 0.2 micron filter and backcheck valve  (15 drops/mL) Universal spike, roller clamp, Supor® 0.2 micron air eliminating filter, backcheck valve approx. 28 in. above distal end, standard injection site 6 in. above distal end, SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 20 mL, Length: 119 in. (approx. 302 cm)		
	V9902F	50

NOTE: To utilize the Free Flow Protection Feature of this set, the height of the IV container should not exceed 40 in. (102 cm) above the backcheck valve. Not intended for gravity infusion.