



Bordered Gauze is an absorptive dressing consisting of three layers: low adherent layer protects the wound surface, absorbent gauze layer absorbs exudate, and a non-woven adhesive tape holds the dressing in place and maintains a moist wound environment. Soft and flexible. Easily conforms to difficult to dress wound sites. Gentle on skin. Ideal for daily use as a primary or secondary dressing.

Bordered Gauze

Adhesive Gauze Dressing

DermaRite's Bordered Gauze dressing has an adhesive border of non-woven cloth tape and comes individually wrapped. The absorptive gauze pad has a low-adherent contact surface that is gentle to the wound bed.

Benefits:

- Non-linting low adherent gauze protects wound bed
- Adhesive backing helps prevent external contamination
- Ideal secondary dressing for daily use
- Latex free

Indications:

For the management of acute or chronic partial and full thickness wounds/ulcers with minimal to moderate exudate

Product	Size	Reorder #	Units
Bordered Gauze	4x4	00255	100/Bag
Bordered Gauze	6x6	00256	100/Bag

DermaRite®