

SenSura® Mio 1-piece Pouching System

The standard in ostomy care



Elastic barrier to fit real body shapes.



Full-circle filter to reduce ballooning.



Neutral gray textile for discretion.

1-piece

SenSura® Mio Deep Convex MAXI Drainable Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE Outlet
- Flex shell is designed for stability and flexibility
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



Transparent 11" (28cm) / 655 mL				Units: 10
Code	Cut	Stoma Size	HCPCS*	
16765	Cut-to-fit	3/8 - 7/8" (10-23 mm)	A4389	
16766	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4389	
16767	Cut-to-fit	3/8 - 1 11/16" (10-43 mm)	A4389	
16770	Pre-cut	1 3/16" (21 mm)	A4389	
16771	Pre-cut	1" (25 mm)	A4389	
16772	Pre-cut	1 7/8" (28 mm)	A4389	
16773	Pre-cut	1 1/4" (31 mm)	A4389	

Opaque 11" (28 cm) / 655 mL				Units: 10
Code	Cut	Stoma Size	HCPCS*	
16755	Cut-to-fit	3/8 - 7/8" (10-23 mm)	A4389	
16756	Cut-to-fit	5/8 - 1 5/16" (15-33 mm)	A4389	
16760	Pre-cut	1 3/16" (21 mm)	A4389	
16761	Pre-cut	1" (25 mm)	A4389	
16762	Pre-cut	1 7/8" (28 mm)	A4389	
16757	Pre-cut	1 1/4" (31 mm)	A4389	
16763	Pre-cut	1 3/8" (35 mm)	A4389	
16764	Pre-cut	1 1/2" (38 mm)	A4389	

SenSura® Mio Deep Convex MAXI Drainable Pouch

- Standard wear
- Belt tabs for use with 4237, 4247
- EasiClose™ WIDE Outlet
- Flex shell is designed for stability and flexibility
- **Without filter**
- Neutral gray, water-resistant textile
- Opaque pouches feature peek-a-boo-inspection window for easy application
- **Not made with natural rubber latex**



Transparent 11" (28 cm) / 655 mL				Units: 10
Code	Cut	Stoma Size	HCPCS*	
16778	Cut-to-fit	3/8 - 1 11/16" (10-43 mm)	A4389	

 Indicates convexity

SenSura® Mio Flat MAXI Closed Pouch

- Standard wear
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window allows easy application
- **Not made with natural rubber latex**



Opaque 8 1/4" (21 cm) / 550 mL				Units: 30
Code	Cut	Stoma Size	HCPCS*	
10881	Cut-to-fit	5/8 - 1 3/4" (15-45 mm)	A4416	
10885**	Cut-to-fit	1 9/16 - 2 1/8" (40-55 mm)	A4416	
10883	Pre-cut	1 3/16" (30 mm)	A4416	
10884	Pre-cut	1 3/8" (35 mm)	A4416	
10885**	Pre-cut	1 9/16" (40 mm)	A4416	

**Because of its larger starter hole, 10885 functions as both a pre-cut and cut-to-fit solution.

SenSura® Mio Flat MIDI Closed Pouch

- Standard wear
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window allows easy application
- **Not made with natural rubber latex**



Opaque 7 1/4" (19 cm) / 340 mL				Units: 30
Code	Cut	Stoma Size	HCPCS*	
10871	Cut-to-fit	5/8 - 1 3/4" (15-45 mm)	A4416	
10872	Pre-cut	1" (25 mm)	A4416	
10873	Pre-cut	1 3/16" (30 mm)	A4416	
10874	Pre-cut	1 3/8" (35 mm)	A4416	

SenSura® Mio Standard Wear MINI Closed Pouch

- **Full-circle filter**
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window allows easy application
- Comfortable and discreet stoma cap alternative
- **Not made with natural rubber latex**



Opaque 5 1/4" (13 cm) / 110 mL				Units: 30
Code	Cut	Stoma Size	HCPCS*	
10870	Cut-to-fit	5/8 - 1 3/4" (15-45 mm)	A4416	

*For important information regarding the use of HCPCS codes listed and code definitions see page 238

The Coloplast story began back in 1954. Our company reflects the passion, ambition and commitment of the people who made it happen. Elise Sørensen, a nurse, invented the first disposable, self-adhesive ostomy pouch, because she was determined to help her sister out of isolation. Aage Louis-Hansen added his engineering ability and entrepreneurial drive and through strong commitment and resilience he founded Coloplast.

Still today, people with intimate healthcare needs often live in isolation.

We are fighting to change that.

For additional information,
contact us at:

1-800-533-0464
info-us@coloplast.com
www.coloplast.us

